**Infectious Bovine Rhinotracheitis Virus (BHV1) gB
Antibody Test Kit**

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For veterinary use only.

Name and Intended Use

IDEXX IBR gB X3 is IDEXX's enzyme immunoassay for the detection of BHV1 specific antibodies to Infectious Bovine Rhinotracheitis Virus (IBR) in individual bovine serum, plasma and milk samples using IBR-gB specific monoclonal antibodies.

General Information

Infectious Bovine Rhinotracheitis is a highly contagious, infectious disease that is caused by Bovine Herpesvirus-1 (BHV1). In addition to causing respiratory disease, this virus can cause conjunctivitis, vulvovaginitis, abortions, encephalitis and generalized systemic infections. Although clinical findings may be highly suggestive of IBR, no real pathognomonic signs are restricted to IBR. Therefore, laboratory confirmation is necessary in order to definitely identify BHV1 infection. Confirmation of exposure to BHV1 via natural infection is facilitated by a measurement of antibody in serum or milk. The enzyme-linked immunosorbent assay (ELISA) for the detection of antibodies against BHV1 in cattle has been shown to correlate with the virus neutralization test (VN), although it can be more sensitive.

Descriptions and Principles

IDEXX IBR gB X3 is an enzyme immunoassay designed to detect the presence of antibodies to IBR/IPV in individual bovine serum, plasma and milk samples. In addition, antibody responses induced by vaccines which contain the glycoprotein B (gB) of BHV1 are detected by this ELISA. A microtitration format has been configured by immobilizing IBR viral antigens on the plate. Upon incubation of the test sample in the antigen-coated well, antibody specific to IBR forms a complex with the immobilized viral antigens. After washing away unbound materials from the wells, a gB-specific monoclonal antibody Horseradish Peroxidase conjugate is added which will not bind to the BHV1 antigen when the antigenic determinant, recognized by the monoclonal antibody, has been occupied (blocked) previously by antibodies in the test sample. Next, unbound conjugate is washed away and a substrate solution is added. In the presence of enzyme, substrate is converted into a product which reacts with the chromogen to generate a blue color. Upon addition of stop solution, a yellow color is generated. The absorbance at a single wavelength of 450 nm [A(450)] or a dual wavelength of 450 nm and 650 nm [A(450/650)] is measured using a spectrophotometer. The blocking percentage of the samples is calculated by using the absorbance [A(450)] or [A(450/650)] obtained with the test sample and the negative control containing no specific antibodies.

Reagents**Volume**

1	BHV-1 Antigen Coated Plate	5	30
2	Positive Control	1 x 2.0 mL	1 x 6.5 mL
3	Negative Control	1 x 2.0 mL	1 x 6.5 mL
4	Conjugate	1 x 60 mL	1 x 350 mL
A	TMB Substrate N.12	1 x 60 mL	1 x 400 mL
B	Stop Solution N.3	1 x 60 mL	1 x 400 mL
C	Wash Concentrate (10X)	1 x 480 mL	3 X 480 mL
Other Components: Zip lock bag		1	1

Note: See table at the end of the insert for a description of symbols used on the insert and labels of this kit.

Storage

Store the reagents at 2–8°C. Reagents are stable until expiration date, provided they have been stored properly.

Materials Required but Not Provided

- Precision micropipettes or multi-dispensing micropipettes
- Disposable pipette tips
- Graduated cylinder for wash solution
- 96-well microplate reader (equipped with single 450 nm filter or dual filters at 450 nm and 650 nm)
- Microplate washer (manual, semi-automatic or automatic system)
- Use only distilled or deionized water for preparation of the reagents used in the test
- Vortex or equivalent
- Microplate covers (lid, aluminium foil or adhesive) and/or humid chamber
- Centrifuge (capacity 2000 x g)
- Microplate shaker
- Incubator capable of maintaining a temperature of +37°C (±3°C)

Precautions and Warnings

- Handle all biological material as potentially infectious.
- Wear protective gloves / protective clothing / eye or face protection when handling samples and reagents.
- Refer to the product Material Safety Data Sheet for additional information.
- See the end of this insert for reagent hazard and precaution warnings.

Laboratory Practices

- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing, and washing throughout this procedure are necessary to maintain precision and accuracy. Use a separate pipette tip for each sample and control.
- Do not expose TMB solution to strong light or any oxidizing agents. Handle TMB solution with clean glass or plastic ware.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional, and national regulations.
- Care should be taken to prevent contamination of kit components. Do not pour unused reagents back into containers.
- Do not use kit past expiration date.

Preparation of Wash Solution

The Wash Concentrate (10X) should be brought to 18–26°C and mixed to ensure dissolution of any precipitated salts. The Wash Concentrate must be diluted 1 to 10 with distilled/deionized water before use (e.g., 30 mL of concentrate plus 270 mL of water per plate to be assayed). When prepared under sterile conditions, the wash solution can be stored for one week at 2–8°C.

Preparation of Samples

Fresh, refrigerated or previously frozen serum or plasma can be tested. Whole milk samples can be used after centrifugation for 15 minutes at 2000 xg or left to stay overnight if refrigerated (2–8°C). No pretreatment is needed for defatted milk.

Test Procedure

All reagents must be allowed to come to 18–26°C before use. Reagents should be mixed by gentle inverting and swirling.

Serum or Plasma Samples

-
- 1 Obtain coated plate(s) and record the sample position. If using partial plates, remove only those wells sufficient for samples to be tested. Place the remaining wells, along with the desiccant, in the extra zip lock bag provided and return to 2–8°C.
-
- 2 Add 50 μ L of reconstituted Wash Solution to each well.
-
- 3 Dispense 50 μ L of Negative Control (NC) into duplicate wells.
-
- 4 Dispense 50 μ L of Positive Control (PC) into duplicate wells.
-
- 5 Dispense 50 μ L of samples into the remaining wells.
-
- 6 Mix the content of the microwells by gently tapping the plate or use a microplate shaker.
-
- 7 Cover the wells and incubate for 2 hours (± 5 min.) at 37°C ($\pm 3^\circ$ C) or overnight (12–18 hours) at 2–8°C. With either option, plates should be hermetically sealed or covered and incubated in a humid chamber. Go to step 8 below.

Milk Samples

- 1** Obtain coated plate(s) and record the sample position. If using partial plates, remove only those wells sufficient for samples to be tested. Place the remaining wells, along with the desiccant, in the extra zip lock bag provided and return to 2–8°C.
 - 2** Dispense 100 µL of Negative Control (NC) into duplicate wells.
 - 3** Dispense 100 µL of Positive Control (PC) into duplicate wells.
 - 4** Dispense 100 µL of defatted milk samples (from underneath the fat layer) into single or duplicate wells.
 - 5** Cover the plate and incubate overnight (12–18 hours) at 2–8°C. Plates should be hermetically sealed or covered and incubated in a humid chamber. Go to step 8 below.
-

Common procedure for serum, plasma, and milk samples

- 8** Remove the solution and wash each well with approximately 300 µL of Wash Solution 5 times. Avoid plate drying between plate washings and prior to the addition of the next reagent. Tap each plate onto absorbent material after the final wash to remove any residual wash fluid.
- 9** Dispense 100 µL of Conjugate into each well.
- 10** Incubate for 1 hour (±5 min.) at 18–26°C.
- 11** Repeat step 8.
- 12** Dispense 100 µL of TMB Substrate N.12 into each well.
- 13** Incubate for 10 minutes (±1 min.) at 18–26°C away from direct light.
- 14** Dispense 100 µL of Stop Solution N.3 into each well.
- 15** Measure and record the absorbance of the samples and controls at 450 nm or using a dual wavelength of 450 nm and 650 nm.

16 Calculation:

Controls

$$NC\bar{x} = \frac{NC1 A(450) + NC2 A(450)}{2}$$

$$PC\bar{x} = \frac{PC1 A(450) + PC2 A(450)}{2}$$

Validity criteria

$$NC\bar{x} \geq 0.500$$

$$PC\bar{x} \text{ Blocking \%} > 80$$

For invalid assays, technique may be suspect and the assay should be repeated following a thorough review of the package insert.

Samples

$$\text{Blocking \%} = 100 \times \frac{NC\bar{x} A(450) - \text{Sample } A(450)}{NC\bar{x} A(450)}$$

The presence or absence of antibody to IBR-gB is determined by the blocking percentage for each sample.

17 Interpretation:

Negative

Suspect

Positive

$$\text{Blocking \%} < 45$$

$$45 \leq \text{Blocking \%} < 55$$

$$\text{Blocking \%} \geq 55$$

Note: IDEXX has instrument and software systems available which calculate means and % values and provide data summaries.

For technical assistance:

IDEXX USA Tel: +1 800 548 9997 or +1 207 556 4895

IDEXX Europe Tel: +800 727 43399

Contact your IDEXX area manager or distributor or visit our website: idexx.com/contactlpd

IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries

© 2016 IDEXX Laboratories, Inc. All rights reserved.

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WARNING / ATTENTION / ATENCIÓN / ATENÇÃO / ACHTUNG

H316 / P332+P313 / EUH208

Conjugate / Positive control / Negative Control — Causes mild skin irritation. If skin irritation occurs: Get medical advice/attention. Contains Proclin. May produce an allergic reaction.

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H315 / H319 / P280 / P332+P313 / P337+P313

TMB Substrate — Causes skin irritation. Causes serious eye irritation. Wear protective gloves / eye protection / face protection. If skin irritation occurs: Get medical advice / attention. If eye irritation persists: Get medical advice / attention.

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Stop solution — Harmful if swallowed. Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. May cause respiratory irritation. Wear protective gloves / eye protection / face protection. If skin irritation or rash occurs: Get medical advice / attention. If eye irritation persists: Get medical advice / attention. Wash contaminated clothing before reuse.

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EUH208

Wash concentrate — Contains Proclin. May produce an allergic reaction.

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IDEXX Laboratories, Inc.
One IDEXX Drive
Westbrook, Maine 04092
USA

Manufacturer
IDEXX Switzerland AG
Stationsstrasse 12
CH-3097 Liebefeld-Bern
Switzerland

EU-Representative
IDEXX Europe B.V.
P.O. Box 1334
2130 EK Hoofddorp
The Netherlands

idexx.com

IDEXX