Avian Influenza Virus Antibody Test Kit, ELISA

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(MultiS-Screen)

Avian Influenza Virus Antibody Test Kit, ELISA (MultiS-Screen)

For veterinary use only.

Name and Intended Use

IDEXX AI MultiS-Screen Test Kit is IDEXX's enzyme immunoassay for the detection of antibody to avian influenza virus (AI) in chicken, turkey, duck, ostrich and goose serum.

General Information

Domestic and wild avian species are affected by avian influenza viruses. The disease is characterized by a wide range of responses that include virtually no clinical signs to high mortality. Respiratory signs are common, along with drop in egg production, greenish diarrhea, bloodstained nasal and oral discharges, and cyanosis and edema of the head, comb and wattle. Because of the variation and severity of clinical symptoms, serological testing produces significant advantages to the detection of infected birds. Monitoring for exposure of a flock to AI is facilitated by the measurement of antibody to AI in serum.

Descriptions and Principles

This assay is designed to measure the relative level of antibody to AI in chicken, turkey, duck, ostrich and goose serum. The assay is performed on 96-well plates that have been coated with avian influenza viral antigen. Upon incubation of the test sample in the coated wells, AI specific antibody forms a complex with the coated antigen. After washing away unbound material, an anti-AI monoclonal antibody enzyme conjugate is added to the wells. In the absence of AI antibodies in the test sample, the conjugate is free to bind to the AI antigen on the plate. Conversely, if there are antibodies to AI present in the sample, the anti-AI conjugate is blocked from binding to the antigen. Unbound conjugate is washed away and enzyme substrate is added. Subsequent color development is inversely proportional to the amount of anti-AI antibodies in the test sample.

Reagents		Volume
1	AI Antigen Coated Plate	5
2	Positive Control — diluted chicken anti-AI serum; preserved with sodium azide	1 x 2.0 mL
3	Negative Control — diluted chicken serum non-reactive to AI; preserved with sodium azide	1 x 2.0 mL
4	Conjugate — anti-AI: HRPO conjugate; preserved with gentamicin and Kathon	1 x 55 mL
5	Sample Diluent —buffer, preserved with sodium azide	1 x 235 mL
A	TMB Substrate	1 x 60 mL
B	Stop Solution	1 x 60 mL
C	Wash Concentrate (10X) — preserved with gentamicin	1 x 235 mL

Note: See table at the end of the insert for a description of symbols used on the insert and labels of this kit.

Storage

Store the reagents at 2–8°C. Reagents are stable until expiration date, provided they have been stored properly.

Materials Required but Not Provided

- Precision micropipettes or multi-dispensing micropipettes
- Disposable pipette tips
- Graduated cylinder for wash solution
- 96-well microplate reader (equipped with 650 nm filter)
- Microplate washer (manual, semi-automatic or automatic system)
- Use only distilled or deionized water for preparation of the reagents used in the test
- Vortex or equivalent

Precautions and Warnings

- Handle all biological material as potentially infectious. The antigen used in the kit reagents may not be completely inactivated.
- Wear protective gloves / protective clothing / eye or face protection when handling samples and reagents.
- Refer to the product Material Safety Data Sheet for additional information.
- See the end of this insert for reagent hazard and precaution warnings.

Laboratory Practices

- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing, and washing throughout this procedure are necessary to maintain precision and accuracy. Use a separate pipette tip for each sample and control.
- Do not expose TMB solution to strong light or any oxidizing agents. Handle TMB solution with clean glass or plastic ware.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional, and national regulations.
- Care should be taken to prevent contamination of kit components. Do not pour unused reagents back into containers.
- Do not use kit past expiration date.

Preparation of Wash Solution

The 10X Wash Concentrate must be brought to 18–26°C and mixed to ensure dissolution of any precipitated salts. The Wash Concentrate must be diluted 1/10 with distilled/deionized water before use (e.g., 30 mL of concentrate plus 270 mL of water per plate to be assayed).

Preparation of Samples

Dilute test samples tenfold (1/10) with sample diluent prior to being assayed (e.g., by diluting 15 μL of sample with 135 μL of sample diluent). **NOTE: DO NOT DILUTE CONTROLS.** Be sure to change tips for each sample. Samples must be thoroughly mixed prior to dispensing into the coated plate.

Test Procedure

All reagents must be allowed to come to 18–26°C before use. Mix reagents by gentle inverting or swirling.

- 1 Obtain antigen-coated plate(s) and record the sample position.

- 2 Dispense 100 μL of UNDILUTED Negative Control (NC) into duplicate wells.

- 3 Dispense 100 μL of UNDILUTED Positive Control (PC) into duplicate wells.

- 4 Dispense 100 μL of DILUTED sample into appropriate wells.

- 5 Incubate for 60 minutes (± 5 minutes) at 18–26°C.

- 6 Remove the solution and wash each well with approximately 350 μL of Wash Solution 3–5 times. Avoid plate drying between plate washings and prior to the addition of the next reagent. Tap each plate onto absorbent material after the final wash to remove any residual wash fluid.

- 7 Dispense 100 μL of Conjugate into each well.

- 8 Incubate for 30 minutes (± 2 minutes) at 18–26°C.

- 9 Repeat step 6.

- 10 Dispense 100 μL of TMB Substrate into each well.

- 11 Incubate for 15 minutes (± 1 minute) at 18–26°C.

- 12 Dispense 100 μL of Stop Solution into each well.

- 13 Measure and record absorbance values at 650nm, A(650).

14 Calculations:

Controls

$$NC\bar{x} = \frac{NC1\ A(650) + NC2\ A(650)}{2}$$

$$PC\bar{x} = \frac{PC1\ A(650) + PC2\ A(650)}{2}$$

Validity criteria

$$NC\bar{x} \geq 0.600$$

$$PC\bar{x} / NC\bar{x} < 0.50$$

For invalid assays, technique may be suspect and the assay should be repeated following a thorough review of the package insert.

Samples

$$S/N = \frac{\text{Sample A}(650)}{NC\bar{x}}$$

The presence or absence of antibody to AI is determined by the sample to negative (S/N) ratio for each sample.

15 Interpretation:

Negative
$S/N \geq 0.50$

Positive
$S/N < 0.50$

Heat inactivation of samples may adversely affect results and should be avoided.

Note: IDEXX has instrument and software systems available which calculate results and provide data summaries.

For technical assistance:

IDEXX USA Tel: +1 800 548 9997 or +1 207 556 4895

IDEXX Europe Tel: +800 727 43399

Contact your IDEXX area manager or distributor or visit our website: idexx.com/contactIpd

U.S. Vet. License No. 313

Product Code: 5004.20

IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries

© 2015 IDEXX Laboratories, Inc. All rights reserved.

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WARNING / ATTENTION / ATENCIÓN / ATENÇÃO / ACHTUNG

H319 / P280 / P337+P313

Positive control / Negative Control / Sample Diluent — Causes skin irritation. Causes serious eye irritation. Wear eye protection. If eye irritation persists: Get medical advice/attention.

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H316 / P332+P313

Stop solution — Causes mild skin irritation. Causes serious eye irritation. Wear eye protection. If skin irritation occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention.

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IDEXX Europe B.V.