**Avian Leukosis Virus Antigen Test Kit**

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For veterinary use only.

## Name and Intended Use

IDEXX ALV Ag is an enzyme immunoassay from IDEXX for the detection of avian leukosis virus Antigen p27 in chicken serum, cloacal swabs or albumin samples.

## General Information

Avian leukosis viruses (ALV) produce a variety of neoplastic diseases including lymphoid leukosis, erythroblastosis, myelocytomatosis, and others.<sup>1</sup> The major gs antigen, p27, is highly conserved across ALV subgroups (A,B,C,D,E and J) and the detection of this antigen is the basis for the ALV-Ag test. Exogenous ALVs can be transmitted horizontally, from bird to bird by direct or indirect contact, or vertically, from an infected hen to progeny. Viremia in the hen is strongly associated with congenital transmission of the virus through shedding into egg albumen. Endogenous leukosis virus sequences are present in the genome of most normal chicken lines.<sup>2</sup>

## Descriptions and Principles

This assay is designed to detect p27, an antigen common to all subgroups of ALV including endogenous viruses. The recommended sample types are light albumin and cloacal swab samples. **While serum has been validated for use on the ALV-Ag test, it is not a recommended sample for the detection of exogenous virus because of potential interference from endogenous sequences.**<sup>3</sup> A microtitration format has been developed in which anti-p27 antibody is coated onto 96-well plates. Sample p27 forms a complex with the coated antibody. After washing away unbound material, an anti-p27:HRPO conjugate is added, which binds to attached p27 in the well. In the final step of the assay, unbound conjugate is washed away and enzyme substrate is added to the well. Color development may then be related to the amount of p27 present in the test sample. Because of the viscosity of albumin samples, a modified sample handling/wash protocol is described as a second procedure in this insert.

Reagents		Volume
1	Anti-p27 Antibody Coated Plate	5
2	Positive Control — inactivated virus in buffer with protein stabilizers; preserved with sodium azide	1 x 1.9 mL
3	Negative Control — buffer with protein additives non-reactive for p27; preserved with sodium azide	1 x 1.9 mL
4	Conjugate — (rabbit) anti-p27: HRPO Conjugate; preserved with gentamicin and Proclin	1 x 50 mL
5	Sample Diluent — buffer, with protein stabilizers; preserved with sodium azide	1 x 235 mL
A	TMB Substrate	1 x 60 mL
B	Stop Solution	1 x 60 mL
C	Albumin Wash Concentrate (20X) — (for albumin wash protocol); preserved with gentamicin	1 x 235 mL

**Note:** See table at the end of the insert for a description of symbols used on the insert and labels of this kit.

## Storage

Store the reagents at 2–8°C. Reagents are stable until expiration date, provided they have been stored properly.

## Materials Required but Not Provided

- Precision micropipettes and multi-dispensing micropipettes
- Disposable pipette tips
- Graduated cylinder for wash solution
- 96-well microplate reader (equipped with 650 nm filter)
- Microplate washer (manual, semi-automatic or automatic system)
- Use only distilled or deionized water for preparation of the reagents used in the test
- Vortex or equivalent

## Precautions and Warnings

- Handle all biological material as potentially infectious. The antigen used in the kit reagents may not be completely inactivated.
- Wear protective gloves / protective clothing / eye or face protection when handling samples and reagents.
- Refer to the product Material Safety Data Sheet for additional information.
- See the end of this insert for reagent hazard and precaution warnings.

## Laboratory Practices

- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing, and washing throughout this procedure are necessary to maintain precision and accuracy. Use a separate pipette tip for each sample and control.
- Do not expose TMB solution to strong light or any oxidizing agents. Handle TMB solution with clean glass or plastic ware.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional, and national regulations.
- Care should be taken to prevent contamination of kit components. Do not pour unused reagents back into containers.
- Do not use kit past expiration date.

## Preparation of Samples

**Albumin** — collect light albumin and add directly to the plate without prior dilution. Freeze/thaw the sample to help reduce the viscosity.

**Cloacal swabs** — place cloacal swab into 1 mL culture media or sample diluent and freeze. Prior to testing, bring the sample to 18–26°C and allow coarse material to settle. Pipette 100  $\mu$ L of the supernatant directly onto the ELISA plate.

**Serum** — for general detection of p27, the sample is added directly to the well without prior dilution. Testing of serum samples for exogenous derived p27 is not recommended because of interference from endogenous virus.

## Test Procedure (For Samples Other Than Albumin)

All reagents must be allowed to come to 18–26°C before use. Mix reagents by gentle inverting or swirling.

- 1 Obtain antibody-coated plate(s) and record the sample position.

---
- 2 Dispense 100  $\mu\text{L}$  of UNDILUTED Negative Control (NC) into duplicate wells.

---
- 3 Dispense 100  $\mu\text{L}$  of UNDILUTED Positive Control (PC) into duplicate wells.

---
- 4 Dispense 100  $\mu\text{L}$  of sample into appropriate wells. Samples may be tested in duplicate, but a single well is acceptable. None of the samples are DILUTED for testing.

---
- 5 Cover the plates and incubate for 60 minutes ( $\pm$  5 minutes) at 18–26°C.

---
- 6 Remove the solution and wash each well with approximately 350  $\mu\text{L}$  of distilled or deionized water 3–5 times. Avoid plate drying between plate washings and prior to the addition of the next reagent. Tap each plate onto absorbent material after the final wash to remove any residual wash fluid.

---
- 7 Dispense 100  $\mu\text{L}$  of Conjugate into each well.

---
- 8 Cover the plates and incubate for 60 minutes ( $\pm$  5 minutes) at 18–26°C.

---
- 9 Repeat step 6.

---
- 10 Dispense 100  $\mu\text{L}$  of TMB Substrate into each well.

---
- 11 Cover the plates and incubate for 15 minutes ( $\pm$  1 minute) at 18–26°C.

---
- 12 Dispense 100  $\mu\text{L}$  of Stop Solution into each well.

---
- 13 Measure and record absorbance values at 650nm, A(650).

**Albumin Wash Protocol** – Light albumin samples are sometimes difficult to completely wash from the wells with the standard water wash protocol. A Wash Concentrate (20X) is provided with the kit for use with albumin samples.

**Preparation of Wash Solution** – The Wash Concentrate (20X) should be brought to 18–26°C and mixed to ensure dissolution of any precipitated salts. The albumin wash solution is prepared by diluting the Wash Concentrate (20X) 1 to 20 with distilled/deionized water before use (e.g., 20 mL of Wash Concentrate (20X) added to 380 mL of water for 1 plate).

**Albumin Wash Procedure** – The albumin test procedure is the same as described above with the exception of the wash steps (#6, #9). Aspirate control and test albumin sample wells and wash with approximately 350  $\mu$ L of wash solution. Allow wells to soak for 2 minutes; aspirate liquid contents; repeat four more times without the 2–minute soak.

---

**14** Calculations:

**Controls**

$$NC\bar{x} = \frac{NC1 A(650) + NC2 A(650)}{2}$$

$$PC\bar{x} = \frac{PC1 A(650) + PC2 A(650)}{2}$$

**Validity criteria**

$$PC\bar{x} - NC\bar{x} > 0.200$$

$$NC\bar{x} \leq 0.150$$

For invalid assays, technique may be suspect and the assay should be repeated following a thorough review of the package insert.

**Samples**

$$S/P = \frac{\text{Sample Mean} - NC\bar{x}}{PC\bar{x} - NC\bar{x}}$$

The presence or absence of p27 antigen is determined by relating the A(650) value of the sample to the Positive Control mean. The Positive Control has been standardized and represents significant antigen levels (approx. 10 ng/mL). The relative level of antigen in the sample can be determined by calculating the sample to positive (S/P) ratio.

---

**15** Interpretation:

Negative

Positive

$$S/P \leq 0.20$$

$$S/P > 0.20$$

**Note:** IDEXX has instrument and software systems available which calculate results and provide data summaries.

## **Bibliography**

1. Payne LN, Fadly AM. 1997. Leukosis/Sarcoma Group in Diseases of Poultry, 10th Ed. B.W. Calnek ed. Ames, Iowa: Iowa State University Press; 1997:414-466.
2. Crittenden LB 1981. Exogenous and endogenous leukosis virus genes—A review. Avian Pathology. 1981;10:101-112.
3. Payne LN, Gillespie AM, Howes K. 1993. Unsuitability of chicken sera for detection of exogenous ALV by the group-specific antigen ELISA. Veterinary Record. May 1993:555-557.

### **For technical assistance:**

IDEXX USA Tel: +1 800 548 9997 or +1 207 556 4895

IDEXX Europe Tel: +800 727 43399

Contact your IDEXX area manager or distributor or visit our website: [idexx.com/contactlpd](http://idexx.com/contactlpd)

U.S. Vet. License No. 313

Product Code: 5007.00

IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries

© 2016 IDEXX Laboratories, Inc. All rights reserved.

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IDEXX Europe B.V.