Bovine Pregnancy Test Kit

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For veterinary use only.

Name and Intended Use

The IDEXX Bovine Pregnancy Test Kit is an enzyme-linked immunoassay for the detection of early pregnancy-associated glycoproteins (early PAGs) in serum or EDTA plasma of cattle, serum of sheep and goat and EDTA plasma of water buffalo and bison as a marker for pregnancy.

General Information

Accurate and timely detection of pregnancy in dairy cows is an essential component of today's reproductive management programs. Veterinarians and farmers use early detection of non-pregnant (open) cows to enable faster rebreeding and shorten the calving interval, thereby maximizing milk production and revenue for the farm. Traditionally, pregnancy is determined through rectal palpation of the cow or by transrectal ultrasound (TU). The IDEXX Bovine Pregnancy Test Kit provides a laboratory based method for accurate detection of pregnancy as early as 28 days post-breeding and from 21 days post embryo transfer, which can give veterinarians and dairy farmers an important tool for the early identification of pregnant cows. The test can be used ≥ 60 days post calving.

Sheep, goat, water buffalo and bison samples

Accurate and timely detection of pregnancy in sheep, goat, water buffalo and bison is an essential component of today's reproductive management programs. Veterinarians and farmers use early detection of non-pregnant (open) sheep, goat, water buffalo and bison to enable faster rebreeding, thereby maximizing revenue for the farm. Traditionally, pregnancy is determined through rectal palpation of the water buffalo or the bison or by abdominal ultrasound (sheep, goat). The IDEXX Bovine Pregnancy Test Kit provides a laboratory based method for accurate detection of pregnancy from 28 days post-breeding in goats, from 35 days in sheep, and from 30 days in water buffalo.* For bison no earliest day post breeding is validated. The test can be used ≥ 60 days post calving in buffalo.* For sheep and goat no post kidding claim is given.

*Refer to the IDEXX Bovine Pregnancy Test validation report for performance data.

Descriptions and Principles

The IDEXX Bovine Pregnancy Test Kit is an enzyme-linked immunoassay designed to detect the early presence of pregnancy-associated glycoproteins (PAGs) in serum or EDTA plasma as a marker for pregnancy. A microtiter plate format has been configured by coating an anti-PAG antibody onto the plate. After dilution and incubation of the test sample in the coated well, captured PAG is detected with a PAG-specific antibody (detector solution) and horseradish peroxidase conjugate (HRPO conjugate). Unbound conjugate is washed away and TMB substrate is added to the wells. Color development is proportional to the amount of PAG in the sample.

Reagents		Volume
1	Anti-PAG-Antibody Coated Plate	5
2	Positive Control	1 x 3.5 mL
3	Negative Control	1 x 3.5 mL
4	Conjugate	1 x 60 mL
5	Sample Diluent	1 x 20 mL
A	TMB Substrate N.12	1 x 60 mL
B	Stop Solution N.3	1 x 60 mL
C	Wash Concentrate (10X)	1 x 480 mL
D	Detector Solution	1 x 60 mL
Other Components: Zip lock bag		1

Note: See table at the end of the insert for a description of symbols used on the insert and labels of this kit.

Storage

Store the reagents at 2–8°C. Reagents are stable until expiration date, provided they have been stored properly.

Materials Required but Not Provided

- Precision micropipettes or multi-dispensing micropipettes
- Disposable pipette tips
- Graduated cylinder for wash solution
- 96-well plate reader equipped with 450 nm filter and reference filter (620 nm – 650 nm)
- Microplate washer (manual, semi-automatic or automatic system)
- Use only distilled or deionized water for preparation of the reagents used in the test
- Microplate covers (lid, aluminium foil or adhesive)
- Vortex or equivalent
- Humid Chamber/ Incubator capable of maintaining a temperature of +37°C (±2°C)

Precautions and Warnings

- Handle all biological material as potentially infectious.
- Wear protective gloves / protective clothing / eye or face protection when handling samples and reagents.
- Refer to the product Material Safety Data Sheet for additional information.
- See the end of this insert for reagent hazard and precaution warnings.

Laboratory Practices

- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing, and washing throughout this procedure are necessary to maintain precision and accuracy. Use a separate pipette tip for each sample and control.
- Do not expose TMB solution to strong light or any oxidizing agents. Handle TMB solution with clean glass or plastic ware.

- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional, and national regulations.
- Care should be taken to prevent contamination of kit components. Do not pour unused reagents back into containers.
- Do not use kit past expiration date.

Preparation of Reagents

Wash Solution

The Wash Concentrate (10X) should be brought to 18–26°C and mixed to ensure dissolution of any precipitated salts. The Wash Concentrate (10X) must be diluted 10-fold (1/10) with distilled/deionized water before use (e.g., 30 mL of concentrate plus 270 mL of water per plate to be assayed).

When prepared under sterile conditions, the Wash Solution can be stored for one week at 2–8°C.

Test Procedure

All reagents must be allowed to come to 18–26°C before use. Reagents should be mixed by gentle inverting or swirling.

- 1 Obtain coated plates and record the sample position. If using partial plates, remove only those wells sufficient for samples to be tested. Place the remaining wells, along with the desiccant, in the extra zip lock bag provided and return to 2–8°C.

- 2 Dispense 25 μL of Sample Diluent into wells which will be used for samples and Controls.

- 3 Dispense 100 μL of Negative Control (NC) into two wells of the assay plate.

- 4 Dispense 100 μL of Positive Control (PC) into two wells of the assay plate.

- 5 Dispense 100 μL of samples into appropriate wells.

- 6 Gently tap plate to mix.

- 7 Cover the wells and incubate 60 minutes (± 5 minutes) at 37°C ($\pm 2^\circ\text{C}$). The plates should be tightly sealed to avoid evaporation.

- 8 Remove the solution and wash each well with approximately 300 μL of Wash Solution 3–5 times. Avoid plate drying between plate washings and prior to the addition of the next reagent. Tap each plate onto absorbent material after the final wash to remove any residual wash fluid.

- 9 Dispense 100 μL of Detector Solution into each well.

- 10 Cover the wells and incubate for 30 minutes (± 2 minutes) at 18–26°C.

- 11 Repeat step 8.

- 12 Dispense 100 μL of Conjugate into each well.

- 13 Cover the wells and incubate for 30 minutes (± 2 minutes) at 18–26°C.

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16 Incubate for 15 minutes (\pm 1 minute) at 18–26°C away from direct light.

17 Dispense 100 μ L of Stop Solution N.3 into each well.

18 Measure and record the A(450 nm)–A(REF) for samples and Controls. The reference wavelength is A(620 nm–650 nm).

19 Calculation:

Controls

$$NC\bar{x} = \frac{NC1 A(450-REF) + NC2 A(450-REF)}{2}$$

$$PC\bar{x} = \frac{PC1 A(450-REF) + PC2 A(450-REF)}{2}$$

Validity criteria

$$PC\bar{x} - NC\bar{x} \geq 0.300$$

$$NC\bar{x} \leq 0.150$$

Samples

$$S - N = \text{Sample A}(450-REF) - NC\bar{x}$$

 20 Interpretation (cattle, sheep, goat, water buffalo and bison):

Open

Pregnant

$$S - N < 0.300$$

$$S - N \geq 0.300$$

Pregnant-Recheck (validated in cattle only): Embryonic death is common in early stages of cow pregnancy. PAGs can circulate in cows/heifers for some time after early embryonic death.

An optional recheck of the animal 7 to 10 days later is recommended if the S – N value is equal to or greater than 0.300 and less than 1.000.

IDEXX encourages producers to work closely with their veterinarian to develop a reproductive management program that is appropriate for their operation.

For invalid assays, technique may be suspect and the assay should be repeated following a thorough review of the package insert. **Note:** IDEXX has instrument and software systems available that calculate results and provide data summaries.

For technical assistance:

IDEXX USA Tel: +1 800 548 9997 or +1 207 556 4895

IDEXX Europe Tel: +800 727 43399

Contact your IDEXX area manager or distributor or visit our website: idexx.com/contact/pd

IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

© 2017 IDEXX Laboratories, Inc. All rights reserved.

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WARNING / ATTENTION / ATENCIÓN / ATENÇÃO / ACHTUNG / ATTENZIONE

H319 / P280 / P337+P313

Positive control – Causes serious eye irritation. Wear eye protection. If eye irritation persists: Get medical advice/attention.

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H315 / H319 / P280 / P332+P313 / P337+P313

TMB Substrate – Causes skin irritation. Causes serious eye irritation. Wear protective gloves/eye protection/face protection. If skin irritation occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention.

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Stop solution – Harmful if swallowed. Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. May cause respiratory irritation. Wear protective gloves/eye protection/face protection. If skin irritation or rash occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention. Wash contaminated clothing before reuse.

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EUH208

Detector Solution – Contains Kathon. May produce an allergic reaction.

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IDEXX Laboratories, Inc.
One IDEXX Drive
Westbrook, Maine 04092
USA

Manufacturer
IDEXX Switzerland AG
Stationsstrasse 12
CH-3097 Liebefeld-Bern
Switzerland

EU-Representative
IDEXX Europe B.V.
P.O. Box 1334
2130 EK Hoofddorp
The Netherlands

idexx.com

IDEXX