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Test With Confidence™

IDEXX

Schmallenberg Virus Antibody Test Kit

For veterinary use only.

Name and Intended Use

IDEXX Schmallenberg Ab Test Kit provides a rapid, simple, sensitive and specific method for detecting antibodies against Schmallenberg virus (SBV) and other Simbu serogroup viruses in individual serum and plasma samples from cattle, sheep and goats.

Descriptions and Principles

Microtiter plates are supplied pre-coated with purified Schmallenberg virus nucleoprotein. Dilutions of the samples to be tested are incubated in the wells of these plates. Any antibody specific for Schmallenberg virus or other Simbu serogroup viruses binds to the antigen in the wells and forms an antigen/antibody complex on the plate well surface. Unbound material is removed from the wells by washing.

Anti-ruminant conjugate is added, which binds to the antibodies complexed with recombinant antigen. Unbound conjugate is removed by washing, and the TMB containing substrate is added to the wells. The degree of colour that develops (optical density measured at 450 nm), is directly proportional to the amount of antibody specific for Schmallenberg virus or other Simbu serogroup viruses present in the sample. The result is obtained by comparing the optical density (OD) of the samples with the OD of the Positive Control.

Reagents

Volume

1	SBV Antigen Coated Plate	2
2	Positive Control	1 x 0.6 mL
3	Negative Control	1 x 0.6 mL
4	Conjugate	1 x 24 mL
5	Sample Diluent	1 x 35 mL
A	TMB Substrate N.12	1 x 20 mL
B	Stop Solution N.3	1 x 20 mL
C	Wash Concentrate (10X)	1 x 125 mL
Other Components: Zip lock bag		1

Note: See table at the end of the insert for a description of symbols used on the insert and labels of this kit.

Storage

Store the reagents at 2–8°C. Reagents are stable until expiration date, provided they have been stored properly.

Materials Required but Not Provided

- Precision micropipettes or multi-dispensing micropipettes
- Disposable pipette tips
- Graduated cylinder for wash solution
- 96-well microplate reader (equipped with 450 nm filter)
- Microplate washer (manual, semi-automatic or automatic system)
- Use only distilled or deionized water for preparation of the reagents used in the test
- Microplate covers (lid, aluminium foil or adhesive)
- Vortex or equivalent

Precautions and Warnings

- Handle all biological material as potentially infectious.
- Wear protective gloves / protective clothing / eye or face protection when handling samples and reagents.
- Refer to the product Material Safety Data Sheet for additional information.
- See the end of this insert for reagent hazard and precaution warnings.

Laboratory Practices

- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing, and washing throughout this procedure are necessary to maintain precision and accuracy. Use a separate pipette tip for each sample and control.
- Do not expose TMB solution to strong light or any oxidizing agents. Handle TMB solution with clean glass or plastic ware.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional, and national regulations.
- Care should be taken to prevent contamination of kit components. Do not pour unused reagents back into containers.
- Do not use kit past expiration date.

Preparation of Reagents

Wash Solution

Dilute the Wash Concentrate (10X) 1:10 with water (1 part concentrate with 9 parts water, e.g., 100 mL Wash Concentrate (10X) + 900 mL distilled water).

Test Procedure

All reagents must be allowed to come to 18–26°C before use. Reagents should be mixed by gentle inverting or swirling.

- 1 Obtain coated microplates and record the position of each sample. If using partial plates, remove only those wells sufficient for samples to be tested. Place the remaining wells, along with the desiccant, in the extra zip lock bag provided and return to 2–8°C.

- 2 Dispense 90 μL of Sample Diluent into each well.

- 3 Dispense 10 μL of Positive Control (PC) in duplicate wells.

- 4 Dispense 10 μL of Negative Control (NC) in duplicate wells.

- 5 Dispense 10 μL of samples in the remaining wells.

- 6 Mix the content of the microwells by gently tapping the plate or use a microplate shaker.

- 7 Cover the plate and incubate for 60 minutes (± 5 min.) at 18–26°C.

- 8 Remove the solution and wash each well with approximately 300 μL of Wash Solution 3 times. Avoid plate drying between plate washings and prior to the addition of the next reagent. Tap each plate onto absorbent material after the final wash to remove any residual wash fluid.

- 9 Dispense 100 μL of Conjugate into each well.

- 10 Cover the plate and incubate it for 60 minutes (± 5 min.) at 18–26°C.

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- 13 Incubate for 10 minutes (± 1 min.) at 18–26°C.

- 14 Dispense 100 μL of Stop Solution N.3 into each well.

- 15 Read the results using a photometer at a wavelength of 450 nm.

Note: Make sure to read the plates within two hours after the addition of the Stop Solution.

16 Calculation:

Controls

$$NC\bar{x} = \frac{NC1 A(450) + NC2 A(450)}{2}$$

$$PC\bar{x} = \frac{PC1 A(450) + PC2 A(450)}{2}$$

Validity criteria

$$NC\bar{x} \leq 0.300$$

$$PC\bar{x} \leq 2.000$$

$$PC\bar{x} - NC\bar{x} \geq 0.300$$

For invalid assays, technique may be suspect and the assay should be repeated following a thorough review of the package insert.

Samples

$$S/P \% = 100 \times \frac{\text{Sample A(450)} - NC\bar{x}}{PC\bar{x} - NC\bar{x}}$$

17 Interpretation:

Negative

Suspect

Positive

$$S/P \% < 30$$

$$30 \leq S/P \% < 40$$

$$S/P \% \geq 40$$

If a sample remains suspect after a second run, a new sample of the same animal should be collected and analyzed again. If the new sample is again suspect, the epidemiological situation should be considered.

Note: IDEXX has instrument and software systems available which calculate results and provide data summaries.

For technical assistance:

IDEXX USA Tel: +1 800 548 9997 or +1 207 556 4895

IDEXX Europe Tel: +800 727 43399

Contact your IDEXX area manager or distributor or visit our website: idexx.com/contactIpd

IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries

© 2015 IDEXX Laboratories, Inc. All rights reserved.

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WARNING / ATTENTION / ATENCIÓN / ATENÇÃO / ACHTUNG

H318 / P280 / P305+P351+P338 / P310

Sample Diluent – Causes serious eye damage. Wear eye protection/face protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/physician.

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H315 / H319 / P280 / P332+P313 / P337+P313

TMB Substrate – Causes skin irritation. Causes serious eye irritation. Wear protective gloves/eye protection/face protection. If skin irritation occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention.

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H302 / H315 / H317 / H319 / H335 / P280 / P333+P313 / P337+P313 / P363

Stop solution – Harmful if swallowed. Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. May cause respiratory irritation. Wear protective gloves/eye protection/face protection. If skin irritation or rash occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention. Wash contaminated clothing before reuse.

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IDEXX Laboratories, Inc.
One IDEXX Drive
Westbrook, Maine 04092
USA

Manufacturer
IDEXX Switzerland AG
Stationsstrasse 12
CH-3097 Liebefeld-Bern
Switzerland

EU-Representative
IDEXX Europe B.V.
P.O. Box 1334
2130 EK Hoofddorp
The Netherlands

idexx.com

IDEXX