**Bovine Viral Diarrhoea Virus (BVDV) Antibody Test Kit**

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Test With Confidence™

**IDEXX**

## Bovine Viral Diarrhoea Virus (BVDV) Antibody Test Kit

For veterinary use only.

### Name and Intended Use

IDEXX BVDV Total Ab is IDEXX's enzyme immunoassay for the detection of BVDV antibodies in individual serum, plasma and milk samples and in bulk milk samples from bovines.

### General Information

Bovine Viral Diarrhoea Virus (BVDV), Border Disease Virus (BDV) and Hog Cholera Virus (Classical Swine Fever Virus, CSFV) are the 3 member viruses of the genus Pestivirus within the family Flaviviridae. BVDV is one of the most important pathogenic viruses in cattle, causing considerable losses in the dairy and beef industries worldwide. The typical symptoms of the BVDV infection are diarrhoea, fever, followed by reduction in the milk production. The immune suppressive effect of BVDV may potentiate infection by other micro-organisms. The virus crosses the placenta in infected, pregnant cows causing reproductive losses due to abortions, still-born calves or calves that die early in life. Some calves that survive are immunotolerant to the virus and these animals excrete large amounts of infectious virus for their entire lives. It is important to identify these 'carrier' animals to break the cycle of infections in herds. The carrier animals often die of 'mucosal disease' in the first 2 years of life. As a consequence of the in-utero infection, BVDV is a frequent contaminant of biological products, such as vaccines and pharmaceuticals. BDV causes similar disease syndromes in sheep, while CSFV causes serious losses in the pig industry since it is highly pathogenic and can cause widespread deaths. However, all pestiviruses may occur in cattle as well.

### Descriptions and Principles

IDEXX BVDV Total Ab is an indirect enzyme immunoassay designed to detect BVDV antibodies (Ab) in serum, plasma and milk samples from bovines. A microtitration format has been configured by immobilizing BVDV antigen on the plates. BVDV antibodies of the sample are bound to the antigen on the plates. After incubation of the test sample in the well, captured BVDV antibodies are detected by anti-bovine horseradish peroxidase conjugate. Next, unbound conjugate is washed away and a substrate/chromogen solution is added. In the presence of enzyme, Substrate is converted into a product which reacts with the chromogen to generate a blue color. Upon addition of stop solution, a yellow color is generated. The absorbance is measured using a spectrophotometer. The sample to positive ratio is calculated by using the absorbance obtained with the test sample and a positive control, corrected for the absorbance of the negative control. Color development indicates the presence of BVDV antibodies in the test sample (positive result).

Reagents		Volume
1	BVDV Antigen Coated Plate	5
2	Positive Control	1 x 1.0 mL
3	Negative Control	1 x 1.0 mL
4	Conjugate	1 x 60 mL
5	Sample Diluent	1 x 60 mL
A	TMB Substrate N.12	1 x 60 mL
B	Stop Solution N.3	1 x 60 mL
C	Wash Concentrate (10X)	1 x 480 mL
<b>Other Components:</b> Zip lock bag		1

**Note:** See table at the end of the insert for a description of symbols used on the insert and labels of this kit.

## Storage

Store the reagents at 2–8°C. Reagents are stable until expiration date, provided they have been stored properly.

## Materials Required but Not Provided

- Precision micropipettes or multi-dispensing micropipettes
- Disposable pipette tips
- Graduated cylinder for wash solution
- 96-well microplate reader (equipped with 450 nm filter or 450 and 650 nm filters)
- Microplate washer (manual, semi-automatic or automatic system)
- Use only distilled or deionized water for preparation of the reagents used in the test
- Microplate covers (lid, aluminium foil or adhesive)
- Centrifuge (capacity 2000 x g)
- Vortex or equivalent
- Microplate shaker

## Precautions and Warnings

- Handle all biological material as potentially infectious.
- Wear protective gloves / protective clothing / eye or face protection when handling samples and reagents.
- Refer to the product Material Safety Data Sheet for additional information.
- See the end of this insert for reagent hazard and precaution warnings.

## Laboratory Practices

- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing, and washing throughout this procedure are necessary to maintain precision and accuracy. Use a separate pipette tip for each sample and control.
- Do not expose TMB solution to strong light or any oxidizing agents. Handle TMB solution with clean glass or plastic ware.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional, and national regulations.
- Care should be taken to prevent contamination of kit components. Do not pour unused reagents back into containers.
- Do not use kit past expiration date.

## Preparation of Wash Solution

The Wash Concentrate (10X) should be brought to 18–26°C and mixed to ensure dissolution of any precipitated salts. The wash concentrate must be diluted 1 to 10 with distilled/deionized water before use (e.g., 30 mL of concentrate plus 270 mL of water per plate to be assayed). When prepared under sterile conditions, the Wash Solution can be stored for one week at 2–8°C.

## Preparation of Samples

Fresh or frozen serum, plasma or milk can be tested. Whole-milk samples may be used after centrifugation for 15 minutes at 2000 x g or left to stay overnight if refrigerated (2–8°C). No pre-treatment is needed for defatted milk.

## Test Procedure

All reagents must be allowed to come to 18–26°C before use. Reagents should be mixed by gentle inverting or swirling.

### Serum or Plasma Samples

- 1 Obtain coated plate(s) and record the sample position. If using partial plates, remove only those wells sufficient for samples to be tested. Place the remaining wells, along with the desiccant, in the extra zip lock bag provided and return to 2–8°C.

---

- 2 Add 100  $\mu$ L of the Sample Diluent to each well.

---

- 3 Add 25  $\mu$ L Negative Control (NC) into duplicate wells.

---

- 4 Add 25  $\mu$ L Positive Control (PC) into duplicate wells.

---

- 5 Add 25  $\mu$ L samples into remaining wells. Continue at Step 6.

### Milk Samples

- 1 Obtain coated plate(s) and record the sample position. If using partial plates, remove only those wells sufficient for samples to be tested. Place the remaining wells, along with the desiccant, in the extra zip lock bag provided and return to 2–8°C.

---

- 2 Add 100  $\mu$ L of the Sample Diluent only into the wells for Negative Control and for Positive Control.

---

- 3 Add 25  $\mu$ L Negative Control (NC) into duplicate wells.

---

- 4 Add 25  $\mu$ L Positive Control (PC) into duplicate wells.

---

- 5 Dispense 100  $\mu$ L UNDILUTED milk samples (from underneath the fat layer) into remaining wells of the plate. Continue at Step 6.

### Common procedure for serum, plasma, and milk samples

- 6 Mix the content of the microwells by gently tapping the plate or use a microplate shaker.

---

- 7 Incubate for 90 minutes ( $\pm$  5 min.) at 18–26°C or overnight (12–18 hours) at 2–8°C. The overnight protocol is recommended for bulk milk samples. With either option, plates should be hermetically sealed or covered and incubated in a humid chamber.

---

- 8 Remove the solution and wash each well with approximately 300  $\mu$ L of Wash Solution 5 times. Avoid plate drying between plate washings and prior to the addition of the next reagent. Tap each plate onto absorbent material after the final wash to remove any residual wash fluid.

---

- 9 Dispense 100  $\mu$ L of Conjugate into each well.

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13 Incubate 10 minutes ( $\pm$  1 min.) at 18–26°C.

---

14 Dispense 100  $\mu$ L of Stop Solution N.3 into each well.

---

15 Measure and record the absorbance of the samples and controls at 450 nm or using a dual wavelength of 450 nm and 650 nm.

---

16 Calculation:

**Controls**

$$NC\bar{x} = \frac{NC1 A(450) + NC2 A(450)}{2}$$

$$PC\bar{x} = \frac{PC1 A(450) + PC2 A(450)}{2}$$

---

**Validity criteria**

$$NC\bar{x} \leq 0.250$$

$$PC\bar{x} - NC\bar{x} \geq 0.150$$

---

For invalid assays, technique may be suspect and the assay should be repeated following a thorough review of the package insert.

---

**Samples**

$$S/P = \frac{\text{Sample A}(450) - NC\bar{x}}{PC\bar{x} - NC\bar{x}}$$

The presence or absence of BVDV antibodies in the sample is determined by S/P ratio for each sample.

---

17 Interpretation:

**Serum, Plasma or Individual Milk samples**

Negative	Suspect	Positive
$S/P < 0.20$	$0.20 \leq S/P < 0.30$	$S/P \geq 0.30$

---

**Bulk Milk samples**

Negative	Positive
$S/P < 0.20$	$S/P \geq 0.20$

---

**Note:** IDEXX has instrument and software systems available that calculate results and provide data summaries.

**For technical assistance:**

IDEXX USA Tel: +1 800 548 9997 or +1 207 556 4895

IDEXX Europe Tel: +800 727 43399

Contact your IDEXX area manager or distributor or visit our website: [idexx.com/contactlpd](http://idexx.com/contactlpd)

IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries

© 2016 IDEXX Laboratories, Inc. All rights reserved.

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[idexx.com/contactlpd](http://idexx.com/contactlpd)

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## WARNING / ATTENTION / ATENCIÓN / ATENÇÃO / ACHTUNG

H319 / P280 / P337+P313 / EUH208

**Conjugate** – Causes serious eye irritation. Wear eye protection / face protection. If eye irritation persists: Get medical advice / attention. Contains Proclin. May produce an allergic reaction.

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H315 / H319 / P280 / P332+P313 / P337+P313

**TMB Substrate** – Causes skin irritation. Causes serious eye irritation. Wear protective gloves / eye protection / face protection. If skin irritation occurs: Get medical advice / attention. If eye irritation persists: Get medical advice / attention.

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H302 / H315 / H317 / H319 / H335 / P280 / P333+P313 / P337+P313 / P363

**Stop solution** – Harmful if swallowed. Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. May cause respiratory irritation. Wear protective gloves / eye protection / face protection. If skin irritation or rash occurs: Get medical advice / attention. If eye irritation persists: Get medical advice / attention. Wash contaminated clothing before reuse.

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EUH208

**Positive control / Negative control / Wash concentrate** – Contains Proclin. May produce an allergic reaction.

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IDEXX Laboratories, Inc.  
One IDEXX Drive  
Westbrook, Maine 04092  
USA

*Manufacturer*  
IDEXX Switzerland AG  
Stationsstrasse 12  
CH-3097 Liebefeld-Bern  
Switzerland

*EU-Representative*  
IDEXX Europe B.V.  
P.O. Box 1334  
2130 EK Hoofddorp  
The Netherlands

[idexx.com](http://idexx.com)

**IDEXX**