Milk Pregnancy Test

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Milk Pregnancy Test

Name and Intended Use

The IDEXX Milk Pregnancy Test is an enzyme-linked immunoassay for the detection of pregnancy-associated glycoproteins (PAGs) in milk as a marker for pregnancy. This test is validated for use in milk samples from cattle, goats, sheep and water buffaloes.

General Information

Accurate and timely detection of pregnancy in production animals is an essential component of today's reproductive management programs. Veterinarians and farmers use detection of non-pregnant (open) animals to make decisions regarding re-breeding and to manage calving intervals, thereby maximizing milk production and revenue for the farm.

Bovine Samples

Traditionally, pregnancy is determined through rectal palpation of the cow or by transrectal ultrasound (TU). The IDEXX Milk Pregnancy Test provides a laboratory based method for the detection of pregnancy from ≥ 28 days post-breeding*, which can give veterinarians and dairy farmers an important tool for the identification of open cows, and management of calving intervals. The test can be used ≥ 60 days post calving.

Goat, sheep and water buffalo samples

Accurate and timely detection of pregnancy in sheep, goat and water buffalo is an essential component of today's reproductive management programs. Veterinarians and farmers use early detection of non-pregnant (open) sheep, goat and water buffalo to enable faster rebreeding, thereby maximizing revenue for the farm. Traditionally, pregnancy is determined through rectal palpation of the water buffalo or by abdominal ultrasound (sheep, goat). The IDEXX Milk Pregnancy Test provides a laboratory based method for accurate detection of pregnancy in goats from ≥ 28 days, in sheep from ≥ 60 days, and in water buffalo from ≥ 29 days post-breeding*, which can give veterinarians and dairy farmers an important tool for the identification and management of open animals. The test can be used ≥ 40 days post calving in water buffaloes and ≥ 80 days post lambing in sheep.*

*Refer to the IDEXX Milk Pregnancy Test validation report for performance data.

Descriptions and Principles

A microtiter plate format has been configured by coating an anti-PAG antibody onto the plate. After incubation of the test sample in the coated well, captured PAG is detected with a PAG-specific antibody (detector solution) and horseradish peroxidase conjugate (HRPO conjugate). Unbound conjugate is washed away and TMB substrate is added to the wells. Color development is proportional to the amount of PAG in the sample.

Reagents**Volume**

1	Anti-PAG-Antibody Coated Plate	5	30
2	Positive Control	1 x 5.0 mL	1 x 14 mL
3	Negative Control	1 x 5.0 mL	1 x 14 mL
4	Conjugate	1 x 60 mL	1 x 350 mL
A	TMB Substrate N.12	1 x 60 mL	1 x 400 mL
B	Stop Solution N.3	1 x 60 mL	1 x 400 mL
C	Wash Concentrate (10X)	1 x 480 mL	3 x 480 mL
D	Detector Solution	1 x 60 mL	1 x 350 mL
Other Components: Zip lock bag		1	1

Note: See table at the end of the insert for a description of symbols used on the insert and labels of this kit.

Storage

Store the reagents at 2–8°C. Reagents are stable until expiration date, provided they have been stored properly.

Materials Required but Not Provided

- Precision micropipettes or multi-dispensing micropipettes
- Disposable pipette tips
- Graduated cylinder for wash solution
- 96-well microplate reader (equipped with single 450 nm filter or dual filters at 620 nm and 650 nm)
- Microplate washer (manual, semi-automatic or automatic system). IDEXX recommends using an automated plate washer and can provide washer settings specific to the Milk Pregnancy Test.
- Use only distilled or deionized water for preparation of the reagents used in the test
- Microplate covers (lid, aluminium foil or adhesive)
- Plate shaker / incubator capable of maintaining a temperature of +37°C ($\pm 2^\circ\text{C}$) and a rotation of 450 rpm (± 50 rpm)

Precautions and Warnings

- Handle all biological material as potentially infectious.
- Wear protective gloves / protective clothing / eye or face protection when handling samples and reagents.
- Refer to the product Material Safety Data Sheet for additional information.
- See the end of this insert for reagent hazard and precaution warnings.

Laboratory Practices

- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing, and washing throughout this procedure are necessary to maintain precision and accuracy. Use a separate pipette tip for each sample and control.
- Do not expose TMB solution to strong light or any oxidizing agents. Handle TMB solution with clean glass or plastic ware.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional, and national regulations.
- Care should be taken to prevent contamination of kit components. Do not pour unused reagents back into containers.
- Do not use kit past expiration date.

Preparation of Reagents

Wash Solution

The Wash Concentrate (10X) should be brought to 18–26°C and mixed to ensure dissolution of any precipitated salts. The Wash Concentrate (10X) must be diluted 10-fold (1/10) with distilled/deionized water before use (e.g., 30 mL of concentrate plus 270 mL of water per plate to be assayed).

When prepared under sterile conditions, the Wash Solution can be stored for one week at 2–8°C.

Preparation of Samples

- Poor quality samples may compromise the accuracy of test results. Bronopol or a similar preservative may be used to maintain sample quality. Prior to testing, milk samples should be checked to ensure that they are not soured or separated and are free from contamination.
- Care should be taken to minimize the likelihood of milk carryover from animal to animal during sample collection, particularly when using samples collected for routine herd recording.
- Preserved or fresh, whole or skim milk can be used. Samples exposed to heat treatment during routine milk component analysis can be used. Samples that have been microwaved cannot be used.

Test Procedure

All reagents must be allowed to come to 18–26°C before use. Reagents should be mixed by gentle inverting or swirling.

- 1 Obtain coated plates and record the sample position. If using partial plates, remove only those wells sufficient for samples to be tested. Place the remaining wells, along with the desiccant, in the extra zip lock bag provided and return to 2–8°C.

- 2 Dispense 150 μL of Negative Control (NC) into two wells of the assay plate.

- 3 Dispense 150 μL of Positive Control (PC) into two wells of the assay plate.

- 4 Dispense 150 μL of Sample into appropriate wells.

- 5 Cover the wells and incubate 120 minutes (± 5 min.) at 37°C ($\pm 2^\circ\text{C}$) with the platform rotating at 450 rpm (± 50 rpm). The plates should be tightly sealed with adhesive cover during this step to avoid evaporation.

- 6 Remove the solution and wash each well with approximately 300 μL of Wash Solution 4 times. Avoid plate drying between plate washings and prior to the addition of the next reagent. Tap each plate onto absorbent material after the final wash to remove any residual wash fluid.

- 7 Dispense 100 μL of Detector Solution into each well.

- 8 Cover the wells and incubate for 30 minutes (± 2 min.) at 18–26°C. The plates should be tightly sealed with adhesive cover during this step to avoid evaporation.

- 9 Repeat step 6.

- 10 Dispense 100 μL of Conjugate into each well.

- 11 Cover the wells and incubate for 30 minutes (± 2 min.) at 18–26°C. The plates should be tightly sealed with adhesive cover during this step to avoid evaporation.

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- 14 Incubate for 20 minutes (± 1 min.) at 18–26°C.

- 15 Dispense 100 μL of Stop Solution N.3 into each well.

- 16 Measure and record the A(450 nm)–A(REF) for samples and Controls. The reference wavelength is A(620 nm–650 nm).

17 Calculation:

Controls

$$NC\bar{x} = \frac{NC1 A(450-REF) + NC2 A(450-REF)}{2}$$

$$PC\bar{x} = \frac{PC1 A(450-REF) + PC2 A(450-REF)}{2}$$

Validity criteria

$$PC\bar{x} - NC\bar{x} \geq 0.500$$

$$NC\bar{x} \leq 0.200$$

For invalid assays, technique may be suspect and the assay should be repeated following a thorough review of the product insert. Pregnant or open (not pregnant) status is determined by the corrected OD values (S – N) for each sample.

Samples

$$S - N = \text{Sample A(450-REF)} - NC\bar{x}$$

18 Interpretation:

Interpretation for bovine (cow) milk samples

S – N Value	Standard Interpretation For Cows ≥ 28 days post breeding	Optional Interpretation [‡] For Cows ≥ 45 days post breeding
< 0.100	Open	Open
≥ 0.100 and < 0.250	Re-check	Pregnant
≥ 0.250	Pregnant	

[‡]Embryonic death is common in early stages of pregnancy. PAGs can circulate in cows for some time after early embryonic death. If the S – N value is greater than or equal to 0.100 45+ days post-breeding, these animals can be considered pregnant. For performance in cows ≥ 45 days post breeding using optional interpretation without a re-check zone, please see the IDEXX Milk Pregnancy Test Validation Report.

Interpretation for milk samples from goat, sheep and water buffalo

S – N Value	Standard Interpretation Goats ≥ 28 , sheep ≥ 60 days, and water buffalo ≥ 29 days post breeding
< 0.100	Open
≥ 0.100 and < 0.250	Re-check
≥ 0.250	Pregnant

Note: IDEXX has instrument and software systems available that calculate results and provide data summaries.

For technical assistance:

IDEXX USA Tel: +1 800 548 9997 or +1 207 556 4895

IDEXX Europe Tel: +800 727 43399

Contact your IDEXX area manager or distributor or visit our website: idexx.com/contactlpd

IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries

© 2017 IDEXX Laboratories, Inc. All rights reserved.

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WARNING / ATTENTION / ATENCIÓN / ATENÇÃO / ACHTUNG / UWAGA

H319 / P280 / P337+P313

Positive control – Causes serious eye irritation. Wear eye protection. If eye irritation persists: Get medical advice/attention.

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H315 / H319 / P280 / P332+P313 / P337+P313

TMB Substrate – Causes skin irritation. Causes serious eye irritation. Wear protective gloves/eye protection/face protection. If skin irritation occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention.

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Stop solution – Harmful if swallowed. Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. May cause respiratory irritation. Wear protective gloves/eye protection/face protection. If skin irritation or rash occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention. Wash contaminated clothing before reuse.

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EUH208

Detector Solution – Contains Proclin. May produce an allergic reaction.

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IDEXX Laboratories, Inc.
One IDEXX Drive
Westbrook, Maine 04092
USA

Manufacturer
IDEXX Switzerland AG
Stationsstrasse 12
CH-3097 Liebefeld-Bern
Switzerland

EU-Representative
IDEXX Europe B.V.
P.O. Box 1334
2130 EK Hoofddorp
The Netherlands

idexx.com

IDEXX