**Bovine Viral Diarrhoea Virus (BVDV) Antigen
Test Kit/Serum Plus**

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For veterinary use only.

Name and Intended Use

IDEXX BVDV Ag/Serum Plus is IDEXX's enzyme immunoassay for the detection of BVDV antigens in bovine serum, plasma, whole blood and ear notch tissue samples.

General Information

Bovine Viral Diarrhoea Virus (BVDV), Border Disease Virus (BDV) and Hog Cholera Virus (Classical Swine Fever Virus, CSFV) are the three member viruses of the genus Pestivirus within the family Flaviviridae. BVDV is one of the most important pathogenic viruses in cattle, causing considerable losses in the dairy and beef industries worldwide. The typical symptoms of the BVDV infection are diarrhoea, fever, followed by reduction in the milk production. The immune suppressive effect of BVDV may potentiate infection by other microorganisms. The virus crosses the placenta in infected, pregnant cows causing reproductive losses due to abortions, stillborn calves or calves that die early in life. Some calves that survive are immunotolerant to the virus and these animals excrete large amounts of infectious virus for their entire lives. It is important to identify these persistently infected animals to break the cycle of infections in herds. The persistently infected animals often die of "mucosal disease" in the first two years of life. As a consequence of the in-utero infection, BVDV is a frequent contaminant of biological products, such as vaccines and pharmaceuticals. BDV causes similar disease syndromes in sheep, while CSFV cause serious losses in the pig industry since it is highly pathogenic and can cause widespread deaths.

Descriptions and Principles

IDEXX BVDV Ag/Serum Plus is an enzyme immunoassay designed to detect BVDV antigens (Ag) in serum, plasma, whole blood and ear notch tissue in cattle. A microtitration format has been configured by immobilizing specific monoclonal antibodies for BVDV (E^{rns}) on the plates. BVDV Ag of the sample is captured on the plates. After incubation of the test sample in the well, captured BVDV Ag is detected by specific antibodies and a horseradish-peroxidase conjugate. Next, unbound conjugate is washed away and a substrate/chromogen solution is added. In the presence of enzyme, substrate is converted into a product which reacts with the chromogen to generate a blue color. Upon addition of the stop solution, a yellow color is generated. The absorbance is measured using a spectrophotometer. The corrected OD value of the samples is calculated by using the absorbance obtained with the test sample and corrected for the absorbance of the negative control.

| Reagents | | Volume | | |
|---------------------------------------|---------------------------------|------------|------------|------------|
| 1 | Anti-ERns mAb Coated Plate | 2 | 5 | 30 |
| 2 | Positive Control | 1 x 1.6 mL | 1 x 2.0 mL | 1 x 6.5 mL |
| 3 | Negative Control | 1 x 1.6 mL | 1 x 2.0 mL | 1 x 6.5 mL |
| 4 | Conjugate | 1 x 25 mL | 1 x 60 mL | 1 x 350 mL |
| 5 | Ear Notch Tissue Soaking Buffer | 1 x 80 mL | 2 x 80 mL | 2 x 480 mL |
| A | TMB Substrate N.12 | 1 x 20 mL | 1 x 60 mL | 1 x 400 mL |
| B | Stop Solution N.3 | 1 x 20 mL | 1 x 60 mL | 1 x 400 mL |
| C | Wash Concentrate (10X) | 1 x 125 mL | 1 x 480 mL | 3 x 480 mL |
| D | Detection solution | 1 x 15 mL | 1 x 30 mL | 1 x 180 mL |
| Other Components: Zip lock bag | | 1 | 1 | — |

Note: See table at the end of the insert for a description of symbols used on the insert and labels of this kit.

Storage

Store the reagents at 2–8°C. Reagents are stable until expiration date, provided they have been stored properly.

Materials Required but Not Provided

- Precision micropipettes or multi-dispensing micropipettes
- Disposable pipette tips
- Graduated cylinder for wash solution
- 96-well microplate reader (equipped with single 450 nm filter or dual filters at 450 nm and 650 nm)
- Microplate washer (manual, semi-automatic or automatic system)
- Use only distilled or deionized water for preparation of the reagents used in the test
- Microplate covers (lid, aluminium foil or adhesive)
- Centrifuge (capacity 2000 x g)
- Vortex or equivalent
- Microplate shaker
- Humid Chamber/Incubator capable of maintaining a temperature of +37°C (±3°C)
- Tubes for soaking of ear notch tissue samples

Precautions and Warnings

- Handle all biological material as potentially infectious.
- Wear protective gloves/protective clothing/eye or face protection when handling samples and reagents.
- Refer to the product Material Safety Data Sheet for additional information.
- See the end of this insert for reagent hazard and precaution warnings.

Laboratory Practices

- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing, and washing throughout this procedure are necessary to maintain precision and accuracy. Use a separate pipette tip for each sample and control.
- Do not expose TMB solution to strong light or any oxidizing agents. Handle TMB solution with clean glass or plastic ware.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional, and national regulations.
- Care should be taken to prevent contamination of kit components. Do not pour unused reagents back into containers.
- Do not use kit past expiration date.

Preparation of Reagents

Wash Solution

The Wash Concentrate (10X) should be brought to 18–26°C and mixed to ensure dissolution of any precipitated salts. The Wash Concentrate (10X) must be diluted 1 to 10 with distilled/deionized water before use (e.g., 30 mL of Wash Concentrate (10X) plus 270 mL of water per plate to be assayed). When prepared under sterile conditions, the Wash Solution can be stored for one week at 2–8°C.

Preparation of Samples

Fresh or frozen serum, plasma, whole blood or ear notch tissue samples can be tested.

Ear Notch Tissue Samples

- Use ear notch tissue plugs (samples) of 2–3 mm diameter in size (e.g., sampled by applying ear tags with attached ear notch tissue sampling device).
- If applicable, ear notch samples can remain in the sampling device for incubation.
Note: Fresh, humid, desiccated or frozen ear notch tissue samples can be tested.
- Add 150–250 μ L of IDEXX Ear Notch Soaking Buffer to the ear notch tissue. Make sure the tissue sample is completely covered with the solution (gently tap or mix).
- Allow the samples to soak in soaking buffer between 12 and 24 hours at 18–26°C or over weekend up to 72 hours at 18–26°C (or at 2–8°C). In order to avoid evaporation, it is recommended to incubate tubes sealed or in a humid chamber.
- Aspirate 50 μ L of the soaking buffer for testing. **Note:** Remaining soaking buffer can be separated from the ear notch tissue sample and stored frozen (–20°C) for later testing or for retesting.

Test Procedure

All reagents must be allowed to come to 18–26°C before use. Reagents should be mixed by gentle inverting or swirling.

- 1 Obtain coated plates and record the sample position. If using partial plates, remove only those wells sufficient for samples to be tested. Place the remaining wells, along with the desiccant, in the extra zip lock bag provided and return to 2–8°C.

- 2 Dispense 50 μ L of the Detection antibodies to each well.

- 3 Dispense 50 μ L Negative Control (NC) into the appropriate wells.

- 4 Dispense 50 μ L Positive Control (PC) into the appropriate wells.

- 5 Dispense 50 μ L samples into the remaining wells.

- 6 Mix the content of the microwells by gently tapping the plate or use a microplate shaker.

- 7 Incubate for 2 hours (\pm 5 min.) at 37°C (\pm 3°C) or overnight (12–18 hours) at 2–8°C. With either option, the plates should be tightly sealed or incubated in a humid chamber using plate covers to avoid any evaporation.

- 8 Remove the solution and wash each well with approximately 300 μ L of Wash Solution 5 times. Avoid plate drying between plate washings and prior to the addition of the next reagent. Tap each plate onto absorbent material after the final wash to remove any residual wash fluid. **Important!** Control carefully that no traces of blood are left on the walls or edges of the wells. Additional 2–3 washes can be necessary to remove the blood before proceeding to the next step.

- 9 Dispense 100 μ L of Conjugate into each well.

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- 13 Incubate for 10 minutes (\pm 1 min.) at 18–26°C.

- 14 Dispense 100 μ L of Stop Solution N.3 into each well.

- 15 Measure and record the absorbance of the samples and controls at 450 nm or using a dual wavelength of 450 nm and 650 nm.

16 Calculation:

Controls

$$NC\bar{x} = \frac{NC1 A(450) + NC2 A(450)}{2}$$

$$PC\bar{x} = \frac{PC1 A(450) + PC2 A(450)}{2}$$

Validity criteria

$$NC\bar{x} \leq 0.250$$

$$PC\bar{x} - NC\bar{x} \geq 0.150$$

For invalid assays, technique may be suspect and the assay should be repeated following a thorough review of the package insert. **Note:** IDEXX has instrument and software systems available that calculate means and S-N and provide data summaries.

Samples

$$S - N = \text{Sample A}(450) - NC\bar{x}$$

The presence or absence of BVDV antigen in the sample is determined by the corrected OD value (S-N) for each sample.

17 Interpretation:

Serum, Plasma and Whole Blood Samples

| Negative | Positive |
|--------------------|-----------------|
| $S - N \leq 0.300$ | $S - N > 0.300$ |

Positive results from this assay are valid for calves of any age. Circulating high titers of maternal BVDV antibodies might interfere with the detection of BVDV antigen in serum, plasma and whole blood. Detection of BVDV antigen in serum, plasma and whole blood samples can be less sensitive after antibody intake via colostrum. "False-negative" results can occur after colostrum intake ("diagnostic gap"). In order to exclude influence of colostrum antibodies, it is recommended to test calves before colostrum intake. Negative results of calves after colostrum intake should be confirmed by re-testing at age of more than 30 days. Refer to regulation in your country if different from this description.

Ear Notch Tissue Samples

| Negative | Suspect | Positive |
|--------------------|----------------------------|-----------------|
| $S - N \leq 0.200$ | $0.200 < S - N \leq 0.300$ | $S - N > 0.300$ |

Suspect samples should be retested by using another 50 µL of the same soaking buffer. If there is less than 50 µL remaining, retesting can be done by soaking the ear notch tissue sample again (see Preparation of Samples). If the sample tests as suspect again, a blood sample should be taken and tested by ELISA with the IDEXX BVDV Ag/Serum Plus, virus isolation or PCR for BVDV. If there is any doubt about the status of available, positive live animal, retest another sample collected 7 to 14 days after the initial sample was collected to confirm persistent infection. We recommend that positive results obtained with whole blood should be confirmed using serum or plasma from the same animal. When testing ear notch samples, there is no diagnostic gap.

For technical assistance:

IDEXX USA Tel: +1 800 548 9997 or +1 207 556 4895

IDEXX Europe Tel: +800 727 43399

Contact your IDEXX area manager or distributor or visit our website: idexx.com/contactlpd

IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

© 2016 IDEXX Laboratories, Inc. All rights reserved.

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EUH208

Positive control / Wash concentrate – Contains Kathon. May produce an allergic reaction.

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H316 / P332+P313 / EUH208

Conjugate – Causes mild skin irritation. If skin irritation occurs: Get medical advice / attention. Contains Kathon. May produce an allergic reaction.

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H319 / P280 / P337+P313 / EUH208

Ear Notch Tissue Soaking Buffer – Causes serious eye irritation. Wear eye protection. If eye irritation persists: Get medical advice / attention. Contains Kathon. May produce an allergic reaction.

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H315 / H319 / P280 / P332+P313 / P337+P313

TMB Substrate – Causes skin irritation. Causes serious eye irritation. Wear protective gloves / eye protection / face protection. If skin irritation occurs: Get medical advice / attention. If eye irritation persists: Get medical advice / attention.

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H316 / H319 / P280 / P332+P313 / P337+P313 / EUH208

Detection Antibodies – Causes mild skin irritation. Causes serious eye irritation. Wear protective gloves / eye protection / face protection. If skin irritation occurs: Get medical advice / attention. If eye irritation persists: Get medical advice / attention. Contains Kathon. May produce an allergic reaction.

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Stop solution – Harmful if swallowed. Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. May cause respiratory irritation. Wear protective gloves/eye protection/face protection. If skin irritation or rash occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention. Take off and wash contaminated clothing before reuse. Dispose of contents in accordance with local/regional/national/international regulations.

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IDEXX Switzerland AG
IDEXX Europe B.V.