**Chicken Anemia Virus Antibody Test Kit**

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Test With Confidence™

**IDEXX**

# Chicken Anemia Virus Antibody Test Kit

For veterinary use only.

## Name and Intended Use

IDEXX CAV is IDEXX's enzyme immunoassay for the detection of antibody to Chicken Anemia Virus (CAV) in chicken serum.

## General Information

Chicken Anemia Virus (CAV) is an important pathogen of poultry and has been found in broilers, breeders and SPF flocks. Virus isolation is difficult and time consuming. Screening for the presence of antibody will indicate exposure to the virus. The enzyme-linked immunosorbent assay has been utilized to detect antibody against CAV. This method is quite useful in large-scale testing of flocks for exposure to CAV. The IDEXX CAV assay can be used as a screen for the presence of Chicken Anemia Virus antibodies or to assess response to vaccination.

## Descriptions and Principles

The IDEXX CAV assay is performed in a microtiter well coated with Chicken Anemia Virus antigen. During the first incubation, CAV antibodies present in the sample react with immobilized antigens. After a wash step, an Anti-CAV monoclonal antibody enzyme conjugate is added to the microwell. If no Anti-CAV antibodies are present in the test sample, the Anti-CAV conjugate is free to react with the CAV antigen. Conversely, if Anti-CAV antibodies are present in the sample, the enzymeconjugated monoclonal antibodies are blocked from reacting with the antigen. Following this incubation period, the unreacted conjugate is removed by washing and a substrate/chromogen solution is added. In the presence of enzyme, the substrate is converted to a product which reacts with the chromophore to generate a blue color. The absorbance at 650nm, A(650), is measured using a spectrophotometer. Results are calculated by dividing the A(650) of the sample by the mean A(650) of the negative control, resulting in an S/N value. The quantity of antibodies to CAV is inversely proportional to the A(650) and thus, to the S/N value. The presence of CAV antibodies indicates previous exposure to Chicken Anemia Virus.

## Reagents

### Volume

	Reagents	Volume
1	CAV Antigen Coated Plate; preserved with Proclin	5
2	Positive Control — diluted chicken anti-CAV serum; preserved with sodium azide	1 x 2.0 mL
3	Negative Control — diluted chicken serum non-reactive to CAV; preserved with sodium azide	1 x 2.0 mL
4	Conjugate — anti-CAV: HRPO conjugate; preserved with gentamicin and Proclin	1 x 50 mL
5	Sample Diluent — preserved with sodium azide	1 x 235 mL
A	TMB Substrate	1 x 60 mL

B	Stop Solution	1 x 60 mL
C	Wash Concentrate (10X) — preserved with gentamicin	1 x 235 mL
<b>Other Components:</b> Zip lock bag.		1

**Note:** See table at the end of the insert for a description of symbols used on the insert and labels of this kit.

## Storage

Store the reagents at 2–8°C. Reagents are stable until expiration date, provided they have been stored properly.

## Materials Required but Not Provided

- Precision micropipettes or multi-dispensing micropipettes
- Disposable pipette tips
- Graduated cylinder for wash solution
- 96-well microplate reader (equipped with 650 nm filter)
- Microplate washer (manual, semi-automatic or automatic system)
- Use only distilled or deionized water for preparation of the reagents used in the test
- Vortex or equivalent

## Precautions and Warnings

- Handle all biological material as potentially infectious. The antigen used in the kit reagents may not be completely inactivated.
- Wear protective gloves / protective clothing / eye or face protection when handling samples and reagents.
- Refer to the product Material Safety Data Sheet for additional information.
- See the end of this insert for reagent hazard and precaution warnings.

## Laboratory Practices

- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing, and washing throughout this procedure are necessary to maintain precision and accuracy. Use a separate pipette tip for each sample and control.
- Do not expose TMB solution to strong light or any oxidizing agents. Handle TMB solution with clean glass or plastic ware.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional, and national regulations.
- Care should be taken to prevent contamination of kit components. Do not pour unused reagents back into containers.
- Do not use kit past expiration date.

## Preparation of Wash Solution

The Wash Concentrate (10X) must be brought to 18–26°C and mixed to ensure dissolution of any precipitated salts. The Wash Concentrate (10X) must be diluted 1/10 with distilled/deionized water before use (e.g. 30 mL of concentrate plus 270 mL of water per plate to be assayed).

## Preparation of Samples

**Note:** Dilution of serum samples may be made with microdilution tubes or directly in the microtiter plate by adding the appropriate amount of serum to the previously pipetted Sample Diluent. Be sure to change tips for each sample. Samples must be thoroughly mixed prior to dispensing into the coated plate. **DO NOT DILUTE THE CONTROLS.**

- A. Field Exposure** – For detection of antibody subsequent to field exposure to Chicken Anemia Virus. Dilute serum samples ten-fold (1/10) with sample diluent prior to being assayed (e.g. diluting 10  $\mu\text{L}$  of sample with 90  $\mu\text{L}$  of sample diluent).
- B. Vaccination** – For the detection of elevated antibody levels (e.g. subsequent to CAV vaccination), dilute serum samples 1/100 (e.g. diluting 10  $\mu\text{L}$  of sample with 990  $\mu\text{L}$  of sample diluent). This dilution will allow adequate differentiation of responses from samples within a vaccinated flock. Responses of vaccination will depend upon dose, route of administration, age of flock, etc.

## Test Procedure

All reagents must be allowed to come to 18–26°C before use. Mix reagents by gentle inverting or swirling.

- 1 Obtain antigen coated plate(s) and record the sample position. If using partial plates, remove only those wells sufficient for samples to be tested. Place the remaining wells, along with the desiccant, in the extra ziplock bag provided and return to 2–8°C.
- 

- 2 Dispense 100  $\mu\text{L}$  of UNDILUTED Negative Control (NC) into duplicate wells.
- 

- 3 Dispense 100  $\mu\text{L}$  of UNDILUTED Positive Control (PC) into duplicate wells.
- 

- 4 Dispense 100  $\mu\text{L}$  of DILUTED sample into appropriate wells. Samples may be tested in duplicate but a single well is acceptable.
- 

- 5 Incubate for 60 minutes ( $\pm$  5 minutes) at 18–26°C.
- 

- 6 Remove the solution and wash each well with approximately 350  $\mu\text{L}$  of Wash Solution 3–5 times. Avoid plate drying between plate washings and prior to the addition of the next reagent. Tap each plate onto absorbent material after the final wash to remove any residual wash fluid.
- 

- 7 Dispense 100  $\mu\text{L}$  of Conjugate into each well.
- 

- 8 Incubate for 30 minutes ( $\pm$  2 minutes) at 18–26°C.
- 

- 9 Repeat step 6.
- 

- 10 Dispense 100  $\mu\text{L}$  of TMB Substrate into each well.
- 

- 11 Incubate for 15 minutes ( $\pm$  1 minute) at 18–26°C.
- 

- 12 Dispense 100  $\mu\text{L}$  of Stop Solution into each well.
- 

- 13 Measure and record absorbance values at 650 nm, A(650).

#### 14 Calculations:

##### Controls

$$NC\bar{x} = \frac{NC1 A(650) + NC2 A(650)}{2}$$

$$PC\bar{x} = \frac{PC1 A(650) + PC2 A(650)}{2}$$

---

##### Validity criteria

$$NC\bar{x} \geq 0.600$$

$$PC\bar{x} / NC\bar{x} \leq 0.50$$

For invalid assays, technique may be suspect and the assay should be repeated following a thorough review of the package insert.

---

##### Samples

$$S/N = \frac{\text{Sample A (650)}}{NC\bar{x}}$$

The presence or absence of antibody to CAV is determined by the sample to negative (S/N) ratio for each sample.

---

#### 15 Interpretation:

##### 1/10 Dilution

Negative

Positive

$$S/N > 0.60$$

$$S/N \leq 0.60$$

---

##### 1/100 Dilution

The immune status of a vaccinated flock is best assessed by monitoring and recording antibody levels in representative samples as a function of time. The resulting flock profiles allow an assessment of the distribution of antibody levels and an analysis of changes in status over time.

**Note:** IDEXX has instrument and software systems available which calculate results and provide data summaries.

#### For technical assistance:

IDEXX USA Tel: +1 800 548 9997 or +1 207 556 4895

IDEXX Europe Tel: +800 727 43399

Contact your IDEXX area manager or distributor or visit our website: [idexx.com/contactlpd](http://idexx.com/contactlpd)

U.S. Vet. License No. 313

Product Code: 50C1.00

IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries

© 2016 IDEXX Laboratories, Inc. All rights reserved.

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IDEXX Europe B.V.