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Test With Confidence™

IDEXX

Infectious Bursal Disease Virus Antibody Test Kit

For veterinary use only.

Name and Intended Use

IDEXX IBD is IDEXX's enzyme immunoassay for the detection of antibody to infectious bursal disease virus (IBD) in chicken serum.

General Information

An assessment of immune status as well as serologic identification of IBD requires a measurement of antibody to IBD in serum. Enzyme immunoassay systems have proven efficacious in the quantification of antibody levels to IBD and facilitate the monitoring of immune status in large flocks.

Descriptions and Principles

This assay is designed to measure the relative level of antibody to IBD in chicken serum. Viral antigen is coated on 96-well plates. Upon incubation of the test sample in the coated well, antibody specific to IBD forms a complex with the coated viral antigens. After washing away unbound material from the wells, a conjugate is added which binds to any attached chicken antibody in the wells. Unbound conjugate is washed away and enzyme substrate is added. Subsequent color development is directly related to the amount of antibody to IBD present in the test sample.

Reagents		Volume
1	IBD Antigen Coated Plate	5
2	Positive Control — diluted chicken anti-IBD serum, preserved with sodium azide	1 x 1.9 mL
3	Negative Control — diluted chicken serum non-reactive to IBD, preserved with sodium azide	1 x 1.9 mL
4	Conjugate — (Goat) anti-Chicken: HRPO Conjugate, preserved with gentamicin and Proclin	1 x 50 mL
5	Sample Diluent — preserved with sodium azide	1 x 235 mL
A	TMB Substrate	1 x 60 mL
B	Stop Solution	1 x 60 mL

Note: See table at the end of the insert for a description of symbols used on the insert and labels of this kit.

Storage

Store the reagents at 2–8°C. Reagents are stable until expiration date, provided they have been stored properly.

Materials Required but Not Provided

- Precision micropipettes and multi-dispensing micropipettes
- Disposable pipette tips
- 96-well microplate reader (equipped with 650 nm filter)
- Microplate washer (manual, semi-automatic or automatic system)
- Use only distilled or deionized water for preparation of the reagents used in the test
- Vortex or equivalent

Precautions and Warnings

- Handle all biological material as potentially infectious. The antigen used in the kit reagents may not be completely inactivated.
- Wear protective gloves / protective clothing / eye or face protection when handling samples and reagents.
- Refer to the product Material Safety Data Sheet for additional information.
- See the end of this insert for reagent hazard and precaution warnings.

Laboratory Practices

- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing, and washing throughout this procedure are necessary to maintain precision and accuracy. Use a separate pipette tip for each sample and control.
- Do not expose TMB solution to strong light or any oxidizing agents. Handle TMB solution with clean glass or plastic ware.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional, and national regulations.
- Care should be taken to prevent contamination of kit components. Do not pour unused reagents back into containers.
- Do not use kit past expiration date.

Preparation of Samples

Dilute test samples five hundred fold (1:500) with sample diluent prior to being assayed (e.g., by diluting 1 μL of sample with 500 μL of Sample Diluent). **NOTE: DO NOT DILUTE CONTROLS.** Be sure to change tips for each sample. Samples must be thoroughly mixed prior to dispensing into the coated plate.

Test Procedure

All reagents must be allowed to come to 18–26°C before use. Mix reagents by gentle inverting or swirling.

- 1 Obtain antigen-coated plate(s) and record the sample position.

- 2 Dispense 100 μL of UNDILUTED Negative Control (NC) into duplicate wells.

- 3 Dispense 100 μL of UNDILUTED Positive Control (PC) into duplicate wells.

- 4 Dispense 100 μL of DILUTED sample into appropriate wells. Samples may be tested in duplicate, but a single well is acceptable.

- 5 Incubate for 30 minutes (\pm 2 minutes) at 18–26°C.

- 6 Remove the solution and wash each well with approximately 350 μL of distilled or deionized water 3–5 times. Avoid plate drying between plate washings and prior to the addition of the next reagent. Tap each plate onto absorbent material after the final wash to remove any residual wash fluid.

- 7 Dispense 100 μL of Conjugate into each well.

- 8 Incubate for 30 minutes (\pm 2 minutes) at 18–26°C.

- 9 Repeat step 6

- 10 Dispense 100 μL of TMB Substrate into each well.

- 11 Incubate for 15 minutes (\pm 1 minutes) at 18–26°C.

- 12 Dispense 100 μL of Stop Solution into each well.

- 13 Measure and record absorbance values at 650nm, A(650).

14 Calculation:

Controls

$$\text{NC}\bar{x} = \frac{\text{NC1 A}(650) + \text{NC2 A}(650)}{2}$$

$$\text{PC}\bar{x} = \frac{\text{PC1 A}(650) + \text{PC2 A}(650)}{2}$$

Validity criteria

$$\text{PC}\bar{x} - \text{NC}\bar{x} > 0.075$$

$$\text{NC}\bar{x} \leq 0.150$$

For invalid assays, technique may be suspect and the assay should be repeated following a thorough review of the package insert.

Samples

$$S/P = \frac{\text{Sample Mean} - NC\bar{x}}{PC\bar{x} - NC\bar{x}}$$

$$\text{Log}_{10} \text{ Titer} = 1.09 (\log_{10} S/P) + 3.36^*$$

*Relates S/P at a 1:500 dilution to an endpoint titer.

The presence or absence of antibody to IBD is determined by relating the A(650) value of the unknown to the Positive Control mean. The Positive Control is standardized and represents significant antibody levels to IBD in chicken serum. The relative level of antibody in the sample is determined by calculating the sample to positive (S/P) ratio. Endpoint titers are calculated using the equation described above.

15 Interpretation:

Negative	Positive
$S/P \leq 0.20$	$S/P > 0.20$

A positive result (titer greater than 396) indicates vaccination or other exposure to IBD. Each laboratory should establish its own criterion for immunity with respect to antibody titer based on correlation of IDEXX IBD to current laboratory test methodologies and on historical antibody responses to specific vaccines and vaccination protocols. The immune status of a flock is best assessed by monitoring and recording antibody titers in representative samples as a function of time. The resulting flock profiles allow an assessment of the distribution of antibody titers and an analysis of changes in titer over time.

Note: IDEXX has instrument and software systems available which calculate results and provide data summaries.

For technical assistance:

IDEXX USA Tel: +1 800 548 9997 or +1 207 556 4895

IDEXX Europe Tel: +800 727 43399

Contact your IDEXX area manager or distributor or visit our website: idexx.com/contactIpd

U.S. Vet. License No. 313

Product Code: 5040.00

IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries

© 2016 IDEXX Laboratories, Inc. All rights reserved.

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IDEXX Europe B.V.