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Test With Confidence™

IDEXX

Influenza A Virus Antibody Test Kit

For veterinary use only.

Name and Intended Use

IDEXX Influenza A Test Kit is an enzyme immunoassay for the detection of antibody to influenza A virus in animal serum.

General Information

Domestic and wild animal species are affected by the influenza A virus. The disease is characterized by a wide range of responses that include virtually no clinical signs to variable mortality. The severity of influenza A will vary with the strain and/or isolate of the virus, the age and the immune status of the animal, the presence of concurrent viral infections and whether influenza A is complicated by secondary infections. Respiratory signs include coughing, nasal and/or ocular discharge, sneezing and dyspnea. Hyperthermia may be observed with associated anorexia, weight loss, lethargy and prostration. Because of the variation and severity of clinical symptoms, serological testing has significant advantages. Monitoring for exposure of a flock or herd to influenza A is facilitated by the measurement of antibody to influenza A in serum.

Descriptions and Principles

This assay is designed to measure the relative level of antibody to influenza A in animal serum. The assay is performed on 96-well plates that have been coated with influenza A viral antigen. Upon incubation of the test sample in the coated wells, influenza A specific antibody forms a complex with the coated antigen. After washing away unbound material, an anti-influenza A nucleoprotein monoclonal antibody enzyme conjugate is added to the wells. In the absence of influenza A antibodies in the test sample, the conjugate is free to bind to the influenza A antigen on the plate. Conversely, if there are antibodies to influenza A present in the sample, the anti-influenza A nucleoprotein conjugate is blocked from binding to the antigen. Unbound conjugate is washed away and enzyme substrate is added. Subsequent color development is inversely proportional to the amount of anti-influenza A antibodies in the test sample.

Reagents

		Volume
1	Influenza A Antigen Coated Plate	5
2	Positive Control	1 x 2.0 mL
3	Negative Control	1 x 2.0 mL
4	Conjugate	1 x 55 mL
5	Dilution Buffer	1 x 235 mL
A	TMB Substrate	1 x 60 mL
B	Stop Solution	1 x 60 mL
C	Wash Concentrate (10X)	1 x 235 mL
Other Components: Zip lock bag		1

Note: See table at the end of the insert for a description of symbols used on the insert and labels of this kit.

Storage

Store the reagents at 2–8°C. Reagents are stable until expiration date, provided they have been stored properly.

Materials Required but Not Provided

- Precision micropipettes or multi-dispensing micropipettes
- Disposable pipette tips
- Graduated cylinder for wash solution
- 96-well microplate reader (equipped with 650 nm filter)
- Microplate washer (manual, semi-automatic or automatic system)
- Use only distilled or deionized water for preparation of the reagents used in the test
- Microplate covers (lid, aluminium foil or adhesive)
- Vortex or equivalent
- Tubes for diluting samples

Precautions and Warnings

- Handle all biological material as potentially infectious.
- Wear protective gloves / protective clothing / eye or face protection when handling samples and reagents.
- Refer to the product Material Safety Data Sheet for additional information.
- See the end of this insert for reagent hazard and precaution warnings.

Laboratory Practices

- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing, and washing throughout this procedure are necessary to maintain precision and accuracy. Use a separate pipette tip for each sample and control.
- Do not expose TMB solution to strong light or any oxidizing agents. Handle TMB solution with clean glass or plastic ware.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional, and national regulations.
- Care should be taken to prevent contamination of kit components. Do not pour unused reagents back into containers.
- Do not use kit past expiration date.

Preparation of Wash Solution

The Wash Concentrate (10X) must be brought to 18–26°C and mixed to ensure dissolution of any precipitated salts. The Wash Concentrate (10X) must be diluted 1/10 with distilled/deionized water before use (e.g., 30 mL of Wash Concentrate (10X) plus 270 mL of water per plate to be assayed). When prepared under sterile conditions, the Wash Solution can be stored for one week at 2–8°C.

Preparation of Samples

Dilute test samples tenfold (1/10) with the Dilution Buffer prior to being assayed (e.g. by diluting 15 μL of sample with 135 μL of Dilution Buffer). **NOTE: DO NOT DILUTE THE CONTROLS.** Samples should be mixed thoroughly before being distributed on the plates.

Test Procedure

All reagents must be allowed to come to 18–26°C before use. Mix reagents by gentle inverting or swirling.

- 1 Obtain coated plate(s) and record the sample position. If using partial plates, remove only those wells sufficient for samples to be tested. Place the remaining wells, along with the desiccant, in the extra zip lock bag provided and return to 2–8°C.

- 2 Dispense 100 μL of UNDILUTED Negative Control (NC) into duplicate wells.

- 3 Dispense 100 μL of UNDILUTED Positive Control (PC) into duplicate wells.

- 4 Dispense 100 μL of DILUTED samples into remaining wells.

- 5 Incubate for 60 minutes (± 5 min.) at 18–26°C.

- 6 Remove the solution and wash each well with approximately 350 μL of Wash Solution 3–5 times. Avoid plate drying between plate washings and prior to the addition of the next reagent. Tap each plate onto absorbent material after the final wash to remove any residual wash fluid.

- 7 Dispense 100 μL of Conjugate into each well.

- 8 Incubate for 30 minutes (± 2 min.) at 18–26°C.

- 9 Repeat step 6.

- 10 Dispense 100 μL of TMB Substrate into each well.

- 11 Incubate 15 minutes (± 1 min.) at 18–26°C.

- 12 Dispense 100 μL of Stop Solution into each well.

- 13 Measure and record the absorbance of the samples and controls at 650 nm.

14 Calculations:

Controls

$$NC\bar{x} = \frac{NC1 A(650) + NC2 A(650)}{2}$$

$$PC\bar{x} = \frac{PC1 A(650) + PC2 A(650)}{2}$$

Validity criteria

$$NC\bar{x} \geq 0.600$$

$$PC\bar{x}/NC\bar{x} < 0.50$$

For invalid assays, technique may be suspect and the assay should be repeated following a thorough review of the package insert.

Samples

$$S/N = \frac{\text{Sample A}(650)}{NC\bar{x}}$$

The presence or absence of antibody to Influenza A is determined by the sample to negative (S/N) ratio for each sample.

15 Interpretation:

Avian Cutoff

Negative
 $S/N \geq 0.50$

Positive
 $S/N < 0.50$

Other Animal Species Cutoff

Negative
 $S/N \geq 0.60$

Positive
 $S/N < 0.60$

Heat inactivation of sera from geese may give a weak false positive test result.

Note: IDEXX has instrument and software systems available which calculate results and provide data summaries.

For technical assistance:

IDEXX USA Tel: +1 800 548 9997 or +1 207 556 4895

IDEXX Europe Tel: +800 727 43399

Contact your IDEXX area manager or distributor or visit our website: idexx.com/contact/tpd

IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries

© 2016 IDEXX Laboratories, Inc. All rights reserved.

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WARNING / ATTENTION / ATENCIÓN / ATENÇÃO / ACHTUNG

H319 / P280 / P337+P313

Positive control / Negative Control / Dilution Buffer – Causes skin irritation. Causes serious eye irritation. Wear eye protection. If eye irritation persists: Get medical advice/attention.

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H316 / H319 / P280 / P332+P313 / P337+P313

Stop solution – Causes mild skin irritation. Causes serious eye irritation. Wear eye protection. If skin irritation occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention.

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IDEXX Montpellier SAS
IDEXX Europe B.V.