Neospora Caninum Antibody Test Kit

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For veterinary use only.

Name and Intended Use

IDEXX *Neospora X2* is IDEXX's enzyme immunoassay for the detection of antibody to *Neospora caninum* in bovine serum.

General Information

Neospora caninum is an apicomplexan protozoal parasite which can cause abortion and neonatal morbidity and mortality in cattle, sheep, goats and horses.¹ *Neospora* has been found worldwide and in many areas is the most common cause of abortions on dairies.^{2,3} Cattle may become infected through the ingestion of oocysts shed in the feces of infected dogs or related canids, the definitive host of the parasite, or by congenital transmission.^{1,4}

An ELISA has been developed for the serological diagnosis of *Neospora* infection in cattle which may be useful in epidemiological studies or in management and control measures for reducing *Neospora* infection within a herd. Since definitive diagnosis is made only by isolation of the organism, the serological status of an animal is only one of many criteria to be considered in overall herd management.

Descriptions and Principles

IDEXX *Neospora X2* is an enzyme immunoassay designed to detect the presence of antibody to *Neospora* in bovine serum. A microtitration format has been developed in which *Neospora* antigens are coated on 96 well plates. Upon incubation of the test sample in the coated well, antibody to *Neospora* forms a complex with the coated antigens. After washing away unbound material from the wells, an anti bovine: horseradish peroxidase conjugate is added which binds to any bovine antibody attached in the wells. In the final step of the assay, unbound conjugate is washed away, and enzyme substrate (hydrogen peroxide) and a chromogen, 3,3', 5,5' tetramethylbenzidine, are added to the wells. Subsequent color development is proportional to the amount of antibody present in the test sample.

Reagents**Volume**

	Reagents	Volume
1	<i>Neospora</i> Antigen Coated Plate	2
2	Positive Control — diluted bovine anti- <i>Neospora</i> serum; preserved with sodium azide	1 x 3,0 mL
3	Negative Control — diluted bovine serum non-reactive to <i>Neospora</i> ; preserved with sodium azide	1 x 3,0 mL
4	Conjugate — anti-bovine IgG: HRPO conjugate; preserved with Proclin	1 x 30 mL
5	Sample Diluent — buffer; preserved with sodium azide	1 x 235 mL

A	TMB Substrate	1 x 60 mL
B	Stop Solution	1 x 60 mL
C	Wash Concentrate (10X) — phosphate/tween wash; preserved with gentamicin	1 x 235 mL
Other Components: Zip lock bag		1

Note: See table at the end of the insert for a description of symbols used on the insert and labels of this kit.

Storage

Store the reagents at 2–8°C. Reagents are stable until expiration date, provided they have been stored properly.

Materials Required but Not Provided

- Precision micropipettes or multi-dispensing micropipettes
- Disposable pipette tips
- Graduated cylinder for wash solution
- 96-well microplate reader (equipped with 620–650 nm filters)
- Microplate washer (manual, semi-automatic or automatic system)
- Use only distilled or deionized water for preparation of the reagents used in the test
- Vortex or equivalent

Precautions and Warnings

- Handle all biological material as potentially infectious. The antigen used in the kit reagents may not be completely inactivated.
- Wear protective gloves / protective clothing / eye or face protection when handling samples and reagents.
- Refer to the product Material Safety Data Sheet for additional information.
- See the end of this insert for reagent hazard and precaution warnings.

Laboratory Practices

- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing, and washing throughout this procedure are necessary to maintain precision and accuracy. Use a separate pipette tip for each sample and control.
- Do not expose TMB solution to strong light or any oxidizing agents. Handle TMB solution with clean glass or plastic ware.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional, and national regulations.
- Care should be taken to prevent contamination of kit components. Do not pour unused reagents back into containers.
- Do not use kit past expiration date.

Preparation of Wash Solution

The Wash Concentrate (10X) must be brought to 18–26°C and mixed to ensure dissolution of any precipitated salts. The Wash Concentrate must be diluted 1/10 with distilled/deionized water before use (e.g. 30 mL of concentrate plus 270 mL of water per plate to be assayed).

Preparation of Samples

Dilute test samples one hundred-fold (1:100) with Sample Diluent (e.g., by diluting 5 μL of sample with 500 μL of Sample Diluent). **DO NOT DILUTE THE CONTROLS.** Be sure to change tips for each sample and record the position of each sample. Samples should be mixed prior to dispensing into the *Neospora* coated plate.

Test Procedure

All reagents must be allowed to come to 18–26°C before use. Mix reagents by gentle inverting or swirling.

- 1 Obtain antigen coated plate(s) and record the sample position. If using partial plates, remove only those wells sufficient for samples to be tested. Place the remaining wells, along with the desiccant, in the extra ziplock bag provided and return to 2–8°C.

- 2 Dispense 100 μL of UNDILUTED Negative Control (NC) into duplicate wells.

- 3 Dispense 100 μL of UNDILUTED Positive Control (PC) into duplicate wells.

- 4 Dispense 100 μL of DILUTED sample into appropriate wells. Samples may be tested in duplicate but a single well is acceptable.

- 5 Incubate for 30 minutes (\pm 2 minutes) at 18–26°C.

- 6 Remove the solution and wash each well with approximately 300 μL of Wash Solution 4 times. Avoid plate drying between plate washings and prior to the addition of the next reagent. Tap each plate onto absorbent material after the final wash to remove any residual wash fluid.

- 7 Dispense 100 μL of Conjugate into each well.

- 8 Incubate for 30 minutes (\pm 2 minutes) at 18–26°C.

- 9 Repeat step 6.

- 10 Dispense 100 μL of TMB Substrate into each well.

- 11 Incubate for 15 minutes (\pm 1 minute) at 18–26°C.

- 12 Dispense 100 μL of Stop Solution into each well.

- 13 Measure and record absorbance values at (620-650 nm).

14 Calculations:

Controls

$$NC\bar{x} = \frac{NC1 A(650) + NC2 A(650)}{2}$$

$$PC\bar{x} = \frac{PC1 A(650) + PC2 A(650)}{2}$$

Validity criteria

$$PC\bar{x} - NC\bar{x} \geq 0.150$$

$$NC\bar{x} \leq 0.200$$

For invalid assays, technique may be suspect and the assay should be repeated following a thorough review of the package insert.

Samples

$$S/P = \frac{\text{Sample Mean} - NC\bar{x}}{PC\bar{x} - NC\bar{x}}$$

The presence or absence of antibody to *Neospora* is determined by sample to positive (S/P) ratio for each sample.

15 Interpretation:

Negative

$$S/P < 0.50$$

Positive

$$S/P \geq 0.50$$

Note: IDEXX has instrument and software systems available which calculate results and provide data summaries.

References

1. Dubey, J.P., Schares G. & Ortega-Mora L.M., 2007. Epidemiology and Control of Neosporosis and *Neospora caninum* Clinical Microbiology Reviews. 20(2): 323-367.
2. Anderson, M.L., et al. 1991. *Neospora*-like Protozoan Infection as a Major Cause of Abortion in California Dairy Cattle. J. Am. Vet. Med. Assoc. 198: 241-244.
3. Paré, J., Thurmond, M.C., Hietala, S.K., 1996. Congenital *Neospora caninum* Infection in Dairy Cattle and Associated Calftlood Mortality. Canadian Journal of Veterinary Research. 60:133-139.
4. Anderson, M. L., et al. 1997. Evidence of vertical transmission of *Neospora* sp. infection in dairy cattle. J. Am. Vet. Med. Assoc. 210:1169-1172.

For technical assistance:

IDEXX USA Tel: +1 800 548 9997 or +1 207 556 4895

IDEXX Europe Tel: +800 727 43399

Contact your IDEXX area manager or distributor or visit our website: idexx.com/contactIpd

U.S. Vet. License No. 313

Product Code: 5N05.00

IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries

© 2016 IDEXX Laboratories, Inc. All rights reserved.

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Références

1. Dubey, J.P., Schares G. & Ortega-Mora L.M., 2007. Epidemiology and Control of Neosporosis and *Neospora caninum* Clinical Microbiology Reviews. 20(2): 323-367.
2. Anderson, M.L., et al. 1991. *Neospora*-like Protozoan Infection as a Major Cause of Abortion in California Dairy Cattle. J. Am. Vet. Med. Assoc. 198: 241-244.
3. Paré, J., Thurmond, M.C., Hietala, S.K., 1996. Congenital *Neospora caninum* Infection in Dairy Cattle and Associated Calving Mortality. Canadian Journal of Veterinary Research. 60:133-139.
4. Anderson, M. L., et al. 1997. Evidence of vertical transmission of *Neospora* sp. infection in dairy cattle. J. Am. Vet. Med. Assoc. 210:1169-1172.

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Referências

1. Dubey, J.P., Schares G. & Ortega-Mora L.M., 2007. Epidemiology and Control of Neosporosis and *Neospora caninum* Clinical Microbiology Reviews. 20(2): 323-367.
2. Anderson, M.L., et al. 1991. *Neospora*-like Protozoan Infection as a Major Cause of Abortion in California Dairy Cattle. J. Am. Vet. Med. Assoc. 198: 241-244.
3. Paré, J., Thurmond, M.C., Hietala, S.K., 1996. Congenital *Neospora caninum* Infection in Dairy Cattle and Associated Calftood Mortality. Canadian Journal of Veterinary Research. 60:133-139.
4. Anderson, M. L., et al. 1997. Evidence of vertical transmission of *Neospora* sp. infection in dairy cattle. J. Am. Vet. Med. Assoc. 210:1169-1172.

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Referencias

1. Dubey, J.P., Schares G. & Ortega-Mora L.M., 2007. Epidemiology and Control of Neosporosis and *Neospora caninum* Clinical Microbiology Reviews. 20(2): 323-367.
2. Anderson, M.L., et al. 1991. *Neospora*-like Protozoan Infection as a Major Cause of Abortion in California Dairy Cattle. J. Am. Vet. Med. Assoc. 198: 241-244.
3. Paré, J., Thurmond, M.C., Hietala, S.K., 1996. Congenital *Neospora caninum* Infection in Dairy Cattle and Associated Calfhood Mortality. Canadian Journal of Veterinary Research. 60:133-139.
4. Anderson, M. L., et al. 1997. Evidence of vertical transmission of *Neospora* sp. infection in dairy cattle. J. Am. Vet. Med. Assoc. 210:1169-1172.

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IDEXX Europe B.V.