Beyond

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Beyond is an extender for the prolonged storage of stallion semen at +5°C or +17°C for 14 or 7 days, respectively. Beyond is a ready to use liquid extender. No further components must be added.

1. Composition

Beyond is a chemically defined synthetic medium that contains antibiotics and a fungicide.

Beyond is a clear to slightly milky solution and without colour, or slightly yellowish. These properties may vary slightly from bottle to bottle.

2. Required Accessories

- We recommend EquiPlus with Amikacin and Penicillin (Ref 13570/0261) as pre-extender for centrifugation. Alternatively, any other milk-based extender with antibiotics may be used.
- optionally CushionFluid, Ref 13580/0001, for centrifugation
- Centrifuge
- Refrigerator adjusted to +5°C
- Transport box for +5°C

For +17°C Storage:

- Semen storage unit/ air conditioned box for semen / refrigerator adjusted to +17°C
- Transport box for +17°C

3. Application

3.1. Preparation

Prepare and warm up the required amount of pre-extender to +30°C to +33°C. Make sure that the pre-extender and the ejaculate are at the same temperature. Warm up **Beyond** to room temperature. No further preparations are required for **Beyond**, it is ready-to-use.

3.2. Recommended application

3.2.1. Preparation of chilled semen

The freshly collected stallion semen is maintained at a temperature of +30°C to +33°C and analysed. After removing the gel, volume, concentration and motility are determined. The filtered gel-free ejaculate is either diluted to 50 million sperm cells/ml or 1:1 with pre-warmed pre-extender. The pre-extender and the semen must have the same temperature at the moment of dilution (+/-1°C difference). The pre-extender is slowly added to the ejaculate and the two parts are mixed gently but thoroughly.

3.2.2. Centrifugation

Centrifugation is performed at 600 g for 15 minutes. After that the supernatant has to be removed completely and the pellet is dissolved in 3-5 ml **Beyond** extender. Ensure, that the pellet and **Beyond** have the same temperature (usually room temperature). To achieve a higher yield of sperm cells, centrifugation may be performed with CushionFluid at up to 1000 g for up to 20 minutes. The extended ejaculate is placed on top of a layer of 3-4 ml of CushionFluid in the centrifugation vial. After centrifugation, the sperm cells are located in a layer between CushionFluid and supernatant. CushionFluid and the supernatant have to be removed.

3.2.3. Final dilution

After centrifugation and resuspension with **Beyond**, the volume, motility and concentration are determined. Alternatively, the concentration is estimated by deducting 20% from the concentration measured in the native ejaculate, to account for the average loss during centrifugation. With these data, the required volume of **Beyond** extender is calculated and then added to the re-suspended pellet (at room temperature) for final dilution.

Calculating the final volume of extended semen:

The optimum dilution rate should be between 50 and 100 million sperm cells/ml after final dilution.

Example:

Data of the freshly collected ejaculate:

Volume of ejaculate = 45 ml
Concentration of semen = 145 million/ml
Total of sperm cells in ejaculate = 6.525 billion (6.525×10^9)
Motility = 60%
→ motile sperm in the ejaculate = 3.915 billion

Data after centrifugation:

Volume = 10 ml
Concentration of semen = 522 million/ml
Total of sperm cells = 5.22 billion
Motility = 60%
→ motile sperm = 3.132 billion

Calculation A (600 million sperm cells per insemination dose):

5 insemination doses can be produced from this ejaculate with 600 million motile semen cells per dose ($3.132 \text{ billion} / 600 \text{ million} = 5.22 \text{ doses}$).
The final volume will be 100 ml with 20 ml volume per insemination dose ($5 \times 20 \text{ ml}$), therefore 90 ml of **Beyond** extender must be added.

Calculation B (50 million sperm cells per ml final dilution):

The re-suspended pellet will be filled up to 104 ml with **Beyond** extender to get a final dilution of 50 million sperm cells/ml ($5.22 \text{ billion sperm cells} / 50 \text{ million sperm cells per ml} = 104.4 \text{ ml}$).

Five doses of 20 ml or 4 doses of 25 ml can be produced from the diluted ejaculate.

3.2.4. Storage of diluted semen

Beyond extender is suitable for the preservation of stallion semen for up to 14 days at +5°C. Semen from stallions known to be very sensitive to cooling, may be stored at +17°C. The increased temperature shortens the time of storage to a maximum of 7 days. Quality control after transport and before insemination is highly recommended.

3.2.5. Insemination

Shortly before insemination, the semen is transferred into a syringe where it should be kept only for a few minutes.

4. Storage, shelf life and important hints

4.1. Storage, shelf life

Beyond is stored at +5°C, the expiration date is printed on the label.

Once opened, store any unused extender at +5°C for up to 3 days. Do not freeze-thaw.

4.2. Important hints

Keep the product out of the reach of children.
The product is not suitable for treatment or consumption by human beings or animals.

Before using **Beyond** extender, make sure that the extender shows the properties described under chapter 1., especially regarding its consistency and appearance. If the extender shows a different consistency or appearance as described, it must not be used.

5. Trouble shooting

A couple of factors are crucial for good semen conservation and insemination success, among others:

- Use of exclusively sperm friendly materials and surfaces for semen and extender (collection gloves, artificial vagina, semen and extender bottles, glassware free of mineral residues).
- Good temperature management, particularly avoiding temperatures below +30°C for undiluted semen and below +5°C for diluted semen, as well as temperature fluctuations.
- Avoiding bacterial contamination (appropriate stallion management and preparations of collection; semen collection in a clean environment, appropriate hygiene protocol in the lab).
- Avoiding sperm toxic factors (water, UV light, dust, detergent).
- Keep the semen in the syringe before insemination as shortly as possible. The usage of syringes for storage or shipping of semen diluted in **Beyond** extender is not recommended.

6. Consulting and training

You can rely on the Minitube scientific experts for advice and technical support concerning all aspects of your semen production and/or artificial insemination process. Workshops and seminars are offered regularly and also on an individual basis. Please refer to our web site for information regarding our current program.

7. Production standard and quality control

7.1. Raw material quality

All raw materials used in Minitube media are produced according to GMP and DIN ISO 9001:2015 norms, and are certified for meeting Ph Eur, BP or USP standards. They are tested and comply with the current international guidelines regulating the quality of pharmaceutical substances.

7.2. Production standard

The mixing and bottling is performed under GMP conditions in a controlled atmosphere room and documented with weighing records. Under the supervision of a team of veterinary specialists, the final product is subjected to a series of chemical and physical tests as well as external spermatological batch tests (in certified laboratories), to ensure that Minitube media will consistently provide semen conservation with high spermatological efficacy over a wide range of climatic and process conditions.

7.3. Certificates

A general quality certificate and the batch certificate are available upon request.

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All rights reserved, including the translation.

Reprint, even in extracts only with a prior written permission. No part of this manual may be reproduced or transmitted in any form or by any means electronic, mechanical, photocopying, without the prior written permission of Minitüb GmbH.

All brand product names are trademarks or registered trademarks of the respective title holder.

Errors and technical alterations excepted.

Use product only in perfect condition and in compliance with the manual. Keep safe for future consulting (according to EN ISO 12100).

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