EquiPlus

Verdünnern für Hengstsamen
Culture Medium for Stallion Semen
Medio de conservación para Semen Equino
Культуральная среда для спермы жеребцов

EquiPlus

13570/0202 (LS, 1 l)
13570/0210 (LS, 100 ml)
13570/0215 (AP, 100 ml)

EquiPlus ohne/without/sin/без AB

13570/0201 (1 l)
13570/0240 (100 ml)

EquiPlus Combi

13570/0221 (G, 100 ml)
13570/0220 (LS, 100 ml)

EquiPlus Combi ohne/without/sin/без AB

13570/0222 (100 ml)

31357/0020

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EquiPlus is a culture medium in which stallion semen can be preserved at approx. +5°C for approx. 3 days. **EquiPlus** is also recommended for centrifugation of stallion semen. **EquiPlus** is supplied in powder form and requires sterile ultrapure water to produce the finished extender. **EquiPlus Combi** comes in a special bottle and contains powder and ultrapure water separately.

1. Composition

EquiPlus is an extender based on skim milk.

The products Ref. 13570/0202, /0210, /0215, /0220 and /0221 are supplemented with antibiotics.

2. Package sizes and alternatives

Ref.	Antibiotics per 1000 ml extender:
13570/0202 13570/0210 13570/0220	0.15 g Lincomycin and 0.25 g Spectinomycin
13570/0215	1 g Amikacin and 0.64 g Penicillin
13570/0221	0.25 g Gentamicin

EquiPlus	1 l	100 ml	without antibiotics	Lincomycin/ Spectinomycin	Amikacin/ Penicillin
13570/0201	X		X		
13570/0202	X			X	
13570/0210		X		X	
13570/0215		X			X
13570/0240		X	X		

EquiPlus Combi	1 l	100 ml	without antibiotics	Gentamicin	Lincomycin/ Spectinomycin
13570/0220		X	X		
13570/0221		X		X	
13570/0222		X			X

3. Application

3.1. Preparation of the final extender

3.1.1. EquiPlus in powder form in sealed bags

The content of an **EquiPlus** bag is a fine, white powder with no particles of other colours and no clumps.

Warm the sterile pure water to a temperature of +30°C to +35°C.

Dissolve the content of one bag of **EquiPlus** in 1000 ml respectively 100 ml purified, sterile water by agitating or stirring.

If EquiPlus without antibiotics is used and an antibiotic is added, the pH of the prepared extender should be checked and adjusted if necessary (e.g. by adding sodium hydrogen carbonate). The pH should be at approx. 6.8.

The ready-to-use extender is liquid, contains no clumps or particles, it has a grey-white colour and it is cloudy. The semen dilution can start after an adaptation time of 15 to 20 minutes.

3.1.2. EquiPlus Combi

EquiPlus Combi has the extender powder included in the lid of the bottle, and the sterile pure water in the bottle itself. To mix the powder with the water, screw and push the upper part of the lid so that the powder is released into the sterile pure water (prewarmed to +30°C to +35°C). Then dissolve the powder completely by gently shaking or agitating the bottle until the extender is homogeneous. The ready-to-use extender is liquid, contains no clumps or particles, it has a grey-white colour and it is cloudy. The semen dilution can start after an adaptation time of 15 to 20 minutes.

3.2. Water quality

The quality of the water used in all the procedures is of a great importance. The water must be purified either by distillation, deionisation, reverse osmosis, or preferably a combination of these methods, and if necessary by additional sterilization. The purified water must have the following characteristics:

Conductivity less than 1 µS/cm

Bacteria level below 1 CFU (colony forming units) per 10 ml

Of utmost importance, all potential sources for bacterial contamination of the prepared water must be identified and eliminated. Unsterile tubing, pipes, or storage containers represent some potential sources for contamination.

If there is no appropriate equipment for water treatment available, bi-distilled and sterile water should be bought from a reliable source.

Never use distilled water for household use!

3.3. Recommended application

3.3.1. Preparation of chilled semen

The freshly collected stallion semen is maintained at a temperature of +30°C to +33°C and analysed. After removing the gel, volume, concentration and motility are determined. If the volume of the ejaculate should be reduced, at first a dilution 1:1 with **EquiPlus** extender and centrifugation must be performed (see 3.3.2). Otherwise the necessary volume of extender is calculated depending on the required number of motile sperm cells per insemination dose. The **EquiPlus** extender and the semen must have the same temperature at the moment of the first dilution (+/-1°C). The extender is slowly added to the filtered, gel-free ejaculate and the 2 fractions are mixed gently but thoroughly in several steps.

After the dilution, the semen is no longer maintained in warm conditions. The recommended cooling rate < 0,3°C per minute, respectively from ambient temperature to +5°C it takes up to 5 hours.

Calculating the final volume of extended semen:

Dilution without centrifugation

The optimum dilution rate should be at least 1:3 (semen:**EquiPlus** extender), or 25 million motile semen cells/ml after final dilution. The recommended motile sperm per insemination dose amounts 600 million. The inseminated volume should be 10-30 ml.

Example:

Data of the freshly collected ejaculate:

Volume of ejaculate = 45 ml

Concentration of semen = 145 million/ml

Total number of sperm = 6.525 billion (6.525×10^9)

Motility = 60%

→ motile sperm in the ejaculate = 3.915 billion

From this ejaculate 6 insemination doses can be obtained, each containing 600 million sperm ($3.915 \text{ billion} / 600 \text{ million} = 6.5 \text{ doses}$). To reach an optimum dilution with 25 million motile sperm/ml, the final volume must be 156 ml ($3.915 \text{ billion} / 25 \text{ million} = 156.6 \text{ ml}$). Each semen dose will amount 26 ml ($156 \text{ ml} / 6 \text{ doses} = 26 \text{ ml}$).

Dilution with centrifugation:

The optimum dilution rate and the number of motile sperm per dose remain equal. The volume of the doses may be reduced (10-15 ml).

Example:

Data of the freshly collected ejaculate:

Volume of ejaculate = 45 ml

Concentration of semen = 145 million/ml

Total number of sperm = 6.525 billion

Motility = 60%

→ motile sperm cells in the ejaculate = 3.915 billion

After centrifugation (see 3.3.2.), the volume is determined, motility evaluated and based on counting with the counting chamber, the new number of sperm cells of the ejaculate is calculated. Alternatively, the new concentration is estimated taking the sperm count of the native ejaculate determined by photometer and adjusted by -20% to account for the loss of cells in centrifugation.

Data after centrifugation:

Volume = 10 ml

Concentration of semen = 522 million/ml

Total number of sperm = 5.22 billion

Motility = 60%

→ motile sperm = 3.132 billion

5 insemination doses can be generated from this ejaculate with 600 million motile semen cells per dose ($3.132 \text{ billion} / 600 \text{ million} = 5.2 \text{ doses}$). The final volume will be 50 ml with 15 ml per insemination dose ($5 \times 15 \text{ ml}$). 65 ml of **EquiPlus** extender must be added to the semen.

EquiPlus is recommended for semen preservation up to 3 days at approx. +5°C. A daily quality analysis until the application is highly recommended.

3.3.2. Centrifugation

The prepared **EquiPlus** extender and the gel-free ejaculate are mixed 1+1. Centrifugation is performed at 500-800 g for 5-10 minutes. After centrifugation the supernatant is removed (a proportion of at least 5% of seminal plasma of the supernatant should remain) and the pellet re-suspended with 3-5 ml of **EquiPlus** extender. To achieve a higher yield of sperm cells it is recommended to work with CushionFluid (Ref. 13580/0001). Therefore, the extended ejaculate is put on top of a layer of 3-4 ml of CushionFluid in the centrifugation vial before centrifugation. Centrifugation is performed at up to 1000 g for up to 20 minutes. After centrifugation, the sperm cells are located in a layer between CushionFluid and supernatant. CushionFluid has to be removed (as well as supernatant) before further semen processing.

4. Storage, shelf life and important hints

4.1. Storage, shelf life

The airproof sealed packages must be stored in cool (below +15°C), dry and dark conditions. The minimum shelf life is printed on each package.

Ready-made extender can be kept up to 3 days if stored in the refrigerator (+2°C - +8°C). It can also be aliquoted and stored at -18°C for approx. 6 months.

4.2. Important hints

**Keep the product out of the reach of children.
Do not touch the powder without gloves and do not breathe in.
The product is not suitable for treatment or consumption by human beings or animals.**

Before each use of the extender, make sure to analyse precisely, whether the extender shows the properties described under 3.1., especially regarding its consistency and aspect. If the extender shows a different consistency or aspect as normal, it must not be used.

4.3. Trouble shooting

A couple of factors are crucial for good semen conservation and insemination success, among others:

- Use of exclusively sperm friendly materials and surfaces for semen and extender (collection gloves, semen and extender bottles, glassware free of mineral residues)
- Good temperature management, particularly avoiding temperatures below +30°C for undiluted semen and below +4°C for diluted semen, as well as temperature fluctuations
- Avoiding bacterial contamination (appropriate stallion management and preparations of collection; semen collection in clean ambient, appropriate hygiene protocol in the lab)
- Avoiding sperm toxic factors (water, UV light, dust, detergent)
- Quality of the purified water (see point 3.2)

4.4. Consulting and training

You can rely on the Minitube scientific experts for advice and technical support concerning all aspects of your semen production and/or artificial insemination process. Workshops and seminars are offered regularly and also on an individual basis. Please refer to our web site for information regarding our current program.

5. Production standard and quality control

5.1. Raw material quality

All raw materials used in Minitube media are produced according to GMP and DIN ISO 9001:2015 norms, and are certified for meeting Ph Eur, BP or USP standards. They are tested and comply with the current international guidelines regulating the quality of pharmaceutical substances.

5.2. Production standard

Production of all Minitube media is performed under GMP conditions in a clean room and documented with weighing records. Under expert veterinary supervision, each batch of the final product is subjected to a series of chemical and physical tests and used under practical conditions to conserve semen. This ensures that Minitube media consistently provide semen conservation with high spermatological efficacy over a wide range of climatic and process conditions.

5.3. Certificates

The general quality certificate and the batch certificate are available upon request.

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All rights reserved, including the translation.
Reprint, even in extracts only with a prior written permission. No part of this manual may be reproduced or transmitted in any form or by any means electronic, mechanical, photocopying, without the prior written permission of Minitüb GmbH.
All brand product names are trademarks or registered trademarks of the respective title holder.
Errors and technical alterations excepted.
Use product only in perfect condition and in compliance with the manual.
Keep safe for future consulting (according to EN ISO 12100).
Translated from German.

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