## IDEXX SNAP® Foal IgG Test Kit

For veterinary use only.

English version

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## SNAP® Foal IgG

The SNAP® Foal IgG Test is an enzyme immunoassay for semi-quantitative detection of immunoglobulin G (IgG) in equine serum, plasma or whole blood. Diluted whole blood, serum or plasma and conjugate reagent are added to the SNAP® device. The device is then activated, releasing reagents stored within the device. Color development on the sample spot of the device is proportional to IgG concentration. Color intensity of the sample spot is compared to color corresponding to sample IgG levels of 400 mg/dL and 800 mg/dL. Comparison of the intensities of color development on the sample spot and the calibrator spots allows an assessment of IgG level in the test sample. Color development on the calibrator spots indicates that the assay reagents are active.

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### Storage

SNAP devices and test reagents are stable until the expiration dates when stored at 2–8°C. **Do not allow SNAP devices or test reagents to freeze.**

### Precautions and Warnings

- Use a separate sample diluent bottle and sample loop for each test. For whole-blood samples, use two separate sample loops for each test. **Do not reuse sample loops.**
- The bioactive spots on the SNAP device are dyed for quality-control purposes. This dye washes out during the test and does not interfere with the test result or interpretation.
- The SNAP device must be in a horizontal position on a flat surface while the test is performed.
- Do not use a SNAP device that has been activated prior to the addition of a sample.
- Do not expose the SNAP device to extreme light after activation.
- Do not mix components from kits with different lot numbers.
- Do not use components past their expiration dates.
- **Conjugate / Sample diluent** – Contains Proclin. May produce an allergic reaction.

### Instructions for Use

- Kit reagents and devices must be at 18–25°C when used.
- Whole blood, serum or plasma can be used.
- Whole blood must be anticoagulated with EDTA or heparin, and can be used either fresh or refrigerated at 2–8°C for up to one week.
- Serum or plasma can be used fresh or stored at 2–8°C for up to 7 days.
- Hemolyzed samples will not affect results. EDTA or heparin in whole-blood and plasma samples will not affect results.

### Test Procedure

**IMPORTANT:** Do not depress the activator until indicated.

1. Remove the cap and dropper tip from the sample diluent bottle.

2. Gently mix whole-blood samples by inverting.

3. Using the plastic sample loops provided, carefully immerse ONLY the loop tip into the sample diluent.

Visually confirm that the loop is filled.

For whole-blood samples, we recommend immersing the loop tip in the sample that remains in the cap of the sample collection container. **Immerse the loop tip only. Do not immerse the loop handle in the sample.**

For serum and plasma, use one loop. **For whole blood, use two separate loops.**

4. Transfer the filled loop by immersing and twirling the loop tip in the bottle of sample diluent.

5. Firmly seat the dropper tip on the sample diluent bottle. Mix thoroughly by inverting five times. Hold the sample diluent bottle vertically and dispose of the first 5–10 drops from the bottle.

6. Place the SNAP device on a flat surface. With the bottle tip ½ to 1 inch directly above the SNAP device, carefully apply one drop of diluted sample directly onto the sample spot in the result window. Visually confirm that the drop of diluted sample has wetted the sample spot completely. If the drop of diluted sample has not wetted the sample spot completely, repeat the sample application using a new SNAP device.

7. Remove the cap from the conjugate bottle and pour its contents into the sample well of the SNAP device. (Some of the contents will remain in the conjugate bottle.) The sample will flow across the result window, reaching the activation circle in 30 to 90 seconds.

(Some conjugate will remain in the sample well at activation.)

**Watch the device carefully for color in the activation circle.**


When color F105T appears in the activation circle, push the activator firmly until it is flush with the device body.

Keep the device horizontal to ensure accurate results.

8. Wait 7 minutes. Visually read the test result.

### Interpreting the Test Result

To determine the test result, read the reaction spots in the result window. To determine the level of IgG in the test sample, compare the color of the sample spot to the color of the two calibrator spots.

### Invalid Test Results

**• Improper Sample Application**  
If the drop of diluted sample applied to the SNAP device does not wet the sample spot completely, nonuniform color development on the sample spot can occur. If this happens, repeat the sample application and test using a new SNAP device.

**• Background**  
If the SNAP device is activated too late (the sample is allowed to flow past the activation circle), high background color may result. If colored background obscures the test result, repeat the test.

**• No Color Development**  
If no color develops on the calibrator spots, repeat the test. If no color develops on the calibrator spots for a second test device, contact IDEXX Technical Support.

**• Color Does Not Appear in the Activation Circle**  
If color does not appear in the activation circle within 90 seconds after the conjugate is added to the SNAP device, repeat the test. If this also happens for a second test device, contact IDEXX Technical Support.

### Symbol descriptions

	Use by date		Authorized representative in the European Community
<b>LOT</b>	Batch code (lot)	<b>EC REP</b>	
<b>SN</b>	Serial number		Consult instructions for use
	Temperature limitation	<b>IVD</b>	In vitro diagnostic
	Manufacturer		Date of manufacture
<b>REF</b>	Catalog number		

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