## GB

Mycoses are infections of fungi. Dermatophytoses are infections in the keratinised skin tissue, i.e. in nails, hair and stratum corneum. Dermatophytoses are caused by one of the following: Microsporum, Trichophyton or Epidermophyton (1).

Dermafyf Test is a media for diagnosing dermatophytoses. In veterinary practice it is mainly Microsporum and Trichophyton that cause dermatophytoses.

Dermafyf Test is based on a colour change from yellow to red when growths of Microsporum or Trichophyton are present. The colour change at normal room temperature (approx. 18-20 °C) will usually take place between two to twelve days. Colour changes after 12 days must be considered unsatisfactory.

Growth of Microsporum or Trichophyton can possibly be verified additionally on colony morphology and the macroconidies in the culture.

### Directions for Use

1. Remove the foil on the Dermafyf Test.
2. Take a small, but representative sample of scale and hair from suspected skin areas and place on the plate.
3. Cover with the plastic lid. (The lid has three ribs which admit sufficient air to the culture.)
4. Record the date and patient reference number.
5. Examine the plate daily for any fungus growth or colour change.

Evaluation of Culture.

It is important to distinguish between dermatophytes and saprophytes. This is done through:

1. **Colony Colour**  
Dermatophytes are white. Saprophytes are brown, grey, green or black.
2. **Colour Change**  
Dermatophytes cause colour change in the test media from yellow to red before or at the same time as the first visible fungus growth is detected. Saprophytes have a well established growth before a possible colour change of the test media takes place.

Colour Change	Colony Colour	Result
None	No growth	Negative
Yellow to red (after 2-12 days)	White	Positive
None	Brown, greyish, green	Negative
Yellow to red (after 12 days)	Brown, greyish, green	Negative

The Dermafyf Test Plate is developed in order to facilitate the transport of the specimen onto the slide before microscopy. Clear tape is pressed carefully against the colony and transferred to a slide on to which f.ex. a drop of Bie's liquid has been placed. During the microscopy the macroconidies (macroaleurospores) can be observed and their type determined. Thus the macroscopic cultivation result is confirmed and at the same time the various Microsporum and Trichophyton types can be determined.

If it is difficult to determine the type of the macroconidies, inoculation of the culture on to a Sabouraud agar plate is recommended as the macroconidies will normally then be very clear (2).

Dermafyf Test plates should be stored in a cool dark place. Do not refrigerate.

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